

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

NEW TESTS			
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.			
Test Code	Test Name	Effective Date	Page #
91932	HPV mRNA E6/E7, Rectal	2/3/2014	2
91768	<i>Streptococcus</i> Group B DNA, PCR with Broth Enrichment	2/3/2014	3
91770	<i>Streptococcus</i> Group B DNA, PCR with Broth Enrichment and Reflex to Susceptibility	2/3/2014	4
92014	<i>Clostridium difficile</i> Culture with Reflex to Toxin	3/10/2014	6
91923	Methylphenidate Metabolite, Quantitative, Urine	3/10/2014	7

TEST CHANGES				
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.				
Test Code	Former Test Code	Test Name	Effective Date	Page #
8369		Hepatitis B Virus DNA, Quantitative, Real Time PCR	3/3/2014	7
822	1345	AST (Aspartate Aminotransferase)	3/10/2014	8
91604		Cardio IQ™ Lipoprotein Fractionation, Ion Mobility	3/10/2014	9
37523		Neutrophil Function, Oxidative Burst	3/10/2014	10
38100	F78G	Casein (f78) IgG	3/17/2014	10
34678	F1G	Egg White (f1) IgG	3/17/2014	11
8927	K82	Latex (k82) IgE	3/17/2014	11
3025	RF266	Mace (f266) IgE	3/17/2014	12
3252	RF301	Persimmon (Kaki Fruit) (f301) IgE	3/17/2014	12
2350	RG203	Salt Grass (g203) IgE	3/17/2014	13
91996	M80	Staphylococcal Enterotoxin A (m80) IgE	3/17/2014	13
91997	M81	Staphylococcal Enterotoxin B (m81) IgE	3/17/2014	14
3412	RM204	<i>Ulocladium chartarum</i> (m204) IgE	3/17/2014	14
30337	F4G	Wheat (f4) IgG	3/17/2014	15
16603	6110	QuantiFERON®-TB Gold, (Draw Site Incubated)	3/31/2014	15

DISCONTINUED TESTS			
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.			
Test Code	Test Name	Effective Date	Page #
8954	<i>Borrelia burgdorferi</i> IgG & IgM Abs [EIA] Plus C6 Peptide	3/10/2014	17
8956	<i>Borrelia burgdorferi</i> IgG & IgM Abs + C6 Peptide w/Rfx IB [CDC]	3/10/2014	17
S52557	Heparin Anti-Xa (Unfractionated Heparin)	3/10/2014	17
8944	Lyme Disease C6 Antibodies, Total, ELISA	3/10/2014	17
S49785NY	Magnesium RBC [623Z][NY]	3/10/2014	17

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

RK82	Allergen - Latex Enhanced IgE	3/17/2014	18
------	-------------------------------	-----------	----

SEND OUTS
Please Note: Not all test codes assigned to each assay are listed in the table of contents.
Please refer to the complete listing on the page numbers indicated.

Test Code	Former Test Code	Test Name	Effective Date	Page #
<u>S47280</u>		Amphetamines Panel, Serum/Plasma	3/3/2014	18
<u>92021</u>		Fluphenazine, Serum/Plasma	3/3/2014	18
<u>S41755</u>		Fluphenazine, Serum/Plasma	3/3/2014	19

New Test Offerings

The following tests will be available through Quest Diagnostics on the dates indicated below.

HPV mRNA E6/E7, Rectal	
Clinical Significance	The presence of E6/E7 messenger RNA from 14 high risk HPV types indicates incorporation of HPV DNA into the host cells. Proteins expressed from E6-E7 polycistronic mRNA alter cellular p53 and retinoblastoma protein functions, leading to disruption of cell-cycle check points and cell genome instability. This information, together with the physician's assessment of cytology history and other risk factors may be used to guide patient management.
Effective Date	2/3/2014
Test Code	91932
CPT Codes	87621
Specimen Requirements	3 mL (1.5 mL minimum) Dacron swab collected in Liquid Cytology (PreservCyt®) Preservative (Thin Prep®)
Reject Criteria	Digene vials; SurePath™ vials
Instructions	To collect an anal-rectal sample, a wetted, non-lubricated Dacron swab is used. The Dacron swab is inserted about 3 cm (or until resistance is met) into the anal canal past the anal verge, into the rectal vault. This is done without visualization of the anal canal. Firm lateral pressure is applied to the swab handle as it is rotated and slowly moved in and out. Slowly withdraw swab from the anal canal. Swish the swab vigorously in PreservCyt® fluid in the ThinPrep® vial. Discard the swab. Cap and tighten the ThinPrep® vial. Follow the same procedure if using a brush to collect an anal-rectal sample. Avoid using cotton swab on a wooden stick because the handle may break and splinter during collection. A swab that is grossly contaminated with feces should be discarded and the collection repeated.
Transport Temperature	Room temperature
Specimen Stability	Room temperature: 30 days Refrigerated: 28 days Frozen: Unacceptable
Set-up/Analytic Time	Set up: Mon, Wed, Fri; Report available 3-6 days
Reference Range	Not detected
Always Message	This test was performed using the APTIMA® HPV Assay (Gen-Probe Inc.). This assay detects E6/E7 viral messenger RNA (mRNA) from 14 high-risk HPV types (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68). The performance characteristics of this assay have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Methodology	Transcription-Mediated Amplification				
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86010331</td> <td>HPV mRNA E6/E7, Rectal</td> </tr> </tbody> </table>	Result Code	Result Name	86010331	HPV mRNA E6/E7, Rectal
Result Code	Result Name				
86010331	HPV mRNA E6/E7, Rectal				

Streptococcus Group B DNA, PCR with Broth Enrichment					
Clinical Significance	Detect Group B <i>Streptococcus</i> (GBS) vaginal/rectal colonization to prevent GBS sepsis and other GBS infections in newborns.				
Effective Date	2/3/2014				
Test Code	91768				
CPT Codes	87081, 87150				
Specimen Requirements	<p>Preferred: Vaginal/Rectal swab in Amies gel (blue-cap) transport medium or Amies liquid (red-cap) transport medium</p> <p>Acceptable: Vaginal/Rectal swab in Liquid Stuart transport medium 2 mL LIM broth tube (pre-incubated)</p>				
Reject Criteria	Specimens in transport media other than those listed; received frozen; LIM broth received at room temperature				
Instructions	<p>In order to obtain an adequate specimen, the procedure for specimen collection must be followed closely. Using recommended swab, vaginal-rectal specimens are collected according to the following procedure:</p> <ol style="list-style-type: none"> 1. Wipe away excessive amount of secretion or discharge from the vaginal area. 2. Carefully insert the swab into the lower one-third part of vagina, and sample secretions from the mucosa. 3. Carefully insert the same swab, approximately 2.5 cm beyond the anal sphincter, and gently rotate to sample anal crypts. 4. Replace the swab in its container. <p>If submitting broth, inoculate LIM broth with vaginal/rectal swab and incubate at 35-37° C for 18-24 hours.</p>				
Transport Temperature	<table border="1"> <tr> <td>Amies transport medium:</td> <td>Room temperature</td> </tr> <tr> <td>LIM broth tube:</td> <td>Refrigerated</td> </tr> </table>	Amies transport medium:	Room temperature	LIM broth tube:	Refrigerated
Amies transport medium:	Room temperature				
LIM broth tube:	Refrigerated				
Specimen Stability	<table border="1"> <tr> <td>Amies transport medium</td> <td>Room temperature and Refrigerated: 48 hours Frozen: Unacceptable</td> </tr> <tr> <td>LIM broth tube</td> <td>Room temperature: Unacceptable Refrigerated: 14 days Frozen: Unacceptable</td> </tr> </table>	Amies transport medium	Room temperature and Refrigerated: 48 hours Frozen: Unacceptable	LIM broth tube	Room temperature: Unacceptable Refrigerated: 14 days Frozen: Unacceptable
Amies transport medium	Room temperature and Refrigerated: 48 hours Frozen: Unacceptable				
LIM broth tube	Room temperature: Unacceptable Refrigerated: 14 days Frozen: Unacceptable				
Set-up/Analytic Time	Set up: Daily; Report available: 2 days				
Reference Range	Not detected				
Always Message	<p>This test was performed using BD GeneOhm™ Streptococcus Group B Real-Time PCR Assay.</p> <p>The performance characteristics of this assay have been determined by Quest Diagnostics Nichols Institute. Performance characteristics refer to the analytical performance of the test.</p>				

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Methodology	Polymerase Chain Reaction, Culture Enhanced		
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano		
CPU Mappings	Result Code	Type	Result Name
	86007404	Prompt-Result	Specimen Source:
	86010442		Group B Streptococcus

Streptococcus Group B DNA, PCR with Broth Enrichment and Reflex to Susceptibility	
Clinical Significance	Detect Group B <i>Streptococcus</i> (GBS) vaginal/rectal colonization to prevent GBS sepsis and other GBS infections in newborns.
Effective Date	2/3/2014
Test Code	91770
CPT Codes	87081, 87150
Specimen Requirements	<p>Preferred: Vaginal/Rectal swab in Amies gel (blue-cap) transport medium or Amies liquid (red-cap) transport medium</p> <p>Acceptable: Vaginal/Rectal swab in Liquid Stuart transport medium is acceptable 2 mL LIM broth tube (pre-incubated)</p>
Reject Criteria	Specimens in transport media other than those listed; received frozen; LIM broth received at room temperature
Instructions	<p>In order to obtain an adequate specimen, the procedure for specimen collection must be followed closely. Using recommended swab, vaginal-rectal specimens are collected according to the following procedure:</p> <ol style="list-style-type: none"> 1. Wipe away excessive amount of secretion or discharge from the vaginal area. 2. Carefully insert the swab into the lower one-third part of vagina, and sample secretions from the mucosa. 3. Carefully insert the same swab, approximately 2.5 cm beyond the anal sphincter, and gently rotate to sample anal crypts. 4. Replace the swab in its container. <p>If submitting broth, inoculate LIM Broth with vaginal/rectal swab and incubate at 35-37° C for 18-24 hours.</p>
Transport Temperature	Amies transport medium: Room temperature
	LIM broth tube: Refrigerated
Specimen Stability	Amies transport medium: Room temperature and Refrigerated: 48 hours Frozen: Unacceptable
	LIM broth tube: Room temperature: Unacceptable Refrigerated: 14 days Frozen: Unacceptable
Set-up/Analytic Time	Set up: Daily; Report available: 3 days
Reference Range	Not detected
Always Message	This test was performed using BD GeneOhm™ Streptococcus Group B Real-Time PCR Assay.

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	The performance characteristics of this assay have been determined by Quest Diagnostics Nichols Institute. Performance characteristics refer to the analytical performance of the test.	
Methodology	Polymerase Chain Reaction, Culture Enhanced	
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano	
CPU Mappings	91770 Streptococcus Group B PCR, Broth Enhanced	
	Result Code	Type
	86007404	Prompt-Result
	86010442	
	<i>This is a true reflex. Please build the unit code below separately. Orderable Reflex: 14653X- Susceptibility, Aerobic Bacteria, MIC</i>	
	Result Code	Type
	24614	Prompt-Result
	26708	Prompt-Result
	26652	
	26653	
	26654	
	26655	
	31019	
	26656	
	26658	
	9991	
	26659	
	26662	
	26660	
	26661	
	26657	
	26663	
	26664	
	26665	
	26666	
	31020	
	31021	
	26667	
	26668	
26645		
26682		
26670		

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	26647	Linezolid
	26704	Meropenem
	26672	Nitrofurantoin
	26675	Oxacillin
	26676	Penicillin
	31022	Piperacillin/Tazobactam
	26679	Tetracycline
	26681	Ticarcillin/Clavulanic ACD
	31023	Tigecycline
	25875	Tobramycin
	26680	Trimethoprim/Sulfamethoxaz
	26649	Vancomycin
	26648	Comments
Additional Information	If GBS by PCR is Detected, an Antibiotic Susceptibility will be added at an additional charge (CPT code(s): 87186 or 87184).	

<i>Clostridium difficile</i> Culture with Reflex to Toxin					
Clinical Significance	<i>C. difficile</i> causes 50-75% of antibiotic associated colitis and greater than 90% of antibiotic associated pseudomembranous colitis. Pathogenic <i>C. difficile</i> produce two potent toxins, toxin A and toxin B, that cause diarrhea and colitis.				
Effective Date	3/10/2014				
Test Code	92014				
CPT Codes	87081				
Specimen Requirements	5 grams or 5 mL (2 grams or 2 mL minimum) unpreserved stool collected in a sterile screw-cap container				
Reject Criteria	Received room temperature; received in Cary-Blair transport; specimen other than liquid or semi-formed stool; stool in preservative or mixed with urine; specimen in wrong transport container.				
Instructions	Toxin is heat labile. Freeze as soon as possible. Submit frozen stool in sterile screw-cap container or use anaerobic transport.				
Transport Temperature	Frozen				
Specimen Stability	Room temperature: Unacceptable Refrigerated: 48 hours Frozen -20°: Not established Frozen -70°: 7 days				
Set-up/Analytic Time	Set up: Daily; Report available: 8-9 days				
Reference Range	<table border="1"> <tr> <td><i>Clostridium difficile</i> Culture</td> <td>Not isolated</td> </tr> <tr> <td><i>Clostridium difficile</i> Toxin B, Qualitative Real-time PCR</td> <td>Not detected</td> </tr> </table>	<i>Clostridium difficile</i> Culture	Not isolated	<i>Clostridium difficile</i> Toxin B, Qualitative Real-time PCR	Not detected
<i>Clostridium difficile</i> Culture	Not isolated				
<i>Clostridium difficile</i> Toxin B, Qualitative Real-time PCR	Not detected				
Always Message	<i>Clostridium difficile</i> Toxin B, Qualitative Real-time PCR:				

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	<p>The performance characteristics of this assay for testing isolates has been determined by Quest Diagnostics Nichols Institute. Performance characteristics refer to the analytical performance of the test.</p> <p>Methodology: BD GeneOhm™ <i>C. difficile</i> Toxin B Gene (tcdB) Real-time PCR Qualitative Assay.</p>														
Methodology	Conventional Culture, Real-Time Polymerase Chain Reaction														
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano														
CPU Mappings	<table border="1"> <tr> <td colspan="2">92014-1 Clostridium difficile Culture</td> </tr> <tr> <td>Result Code</td> <td>Result Name</td> </tr> <tr> <td>85996814</td> <td>C. difficile Culture</td> </tr> <tr> <td colspan="2">*TR 92014-2 Clostridium difficile Toxin B, Qualitative Real-time PCR</td> </tr> <tr> <td>Result Code</td> <td>Result Name</td> </tr> <tr> <td>86010665</td> <td>C. difficile, QL PCR</td> </tr> <tr> <td colspan="2"> <p>*TR (True Reflexing Flag) CPU interface clients: If you are set up to use our True Reflexing option, build the unit code with the TR flag (indicated above) separately.</p> </td> </tr> </table>	92014-1 Clostridium difficile Culture		Result Code	Result Name	85996814	C. difficile Culture	*TR 92014-2 Clostridium difficile Toxin B, Qualitative Real-time PCR		Result Code	Result Name	86010665	C. difficile, QL PCR	<p>*TR (True Reflexing Flag) CPU interface clients: If you are set up to use our True Reflexing option, build the unit code with the TR flag (indicated above) separately.</p>	
92014-1 Clostridium difficile Culture															
Result Code	Result Name														
85996814	C. difficile Culture														
*TR 92014-2 Clostridium difficile Toxin B, Qualitative Real-time PCR															
Result Code	Result Name														
86010665	C. difficile, QL PCR														
<p>*TR (True Reflexing Flag) CPU interface clients: If you are set up to use our True Reflexing option, build the unit code with the TR flag (indicated above) separately.</p>															
Additional Information	If <i>Clostridium difficile</i> Culture is "Isolated" then the <i>Clostridium difficile</i> Toxin B, Qualitative Real-time PCR is performed at an additional charge (CPT code(s): 87493)														

Methylphenidate Metabolite, Quantitative, Urine					
Clinical Significance	The analysis of the methylphenidate metabolite is utilized to monitor compliance with prescribed therapy.				
Effective Date	3/10/2014				
Test Code	91923				
CPT Codes	80299				
Specimen Requirements	20 mL (5 mL minimum) random urine collected in sterile plastic leak-proof container				
Transport Temperature	Room temperature				
Specimen Stability	Room temperature: 5 days Refrigerated: 7 days Frozen: 30 days				
Set-up/Analytic Time	Set-up: Tues-Sat; Report available: 2-3 days				
Reference Range	<100 ng/mL				
Methodology	Liquid Chromatography Tandem Mass Spectrometry				
Assay Category	Laboratory Developed Test				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
CPU Mappings	<table border="1"> <tr> <td>Result Code</td> <td>Result Name</td> </tr> <tr> <td>86010283</td> <td>MethylphenidateMetab,QN,U</td> </tr> </table>	Result Code	Result Name	86010283	MethylphenidateMetab,QN,U
Result Code	Result Name				
86010283	MethylphenidateMetab,QN,U				

Test Changes

The following test changes will be effective on the dates indicated below. **Please note that only the information that is changing appears in this**

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

update. Former test names and test codes have been italicized.

Hepatitis B Virus DNA, Quantitative, Real Time PCR											
Effective Date	3/3/2014										
Test Code	8369										
Reference Range	Hepatitis B Virus DNA: <20 IU/mL Hepatitis B Virus DNA: <1.30 LogIU/mL										
Always Message	<p>This test was performed using the COBAS(R) AmpliPrep/COBAS(R) TaqMan(R) HBV Test, v2.0 (Roche Molecular Systems, Inc.).</p> <p>To convert International Units/mL to copies/mL, use the following conversion factor: 1 IU/mL = 5.82 copies/mL.</p>										
Performing Site	Quest Diagnostics Nichols Institute, Valencia										
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> <th>Unit of Measure</th> </tr> </thead> <tbody> <tr> <td>85997326</td> <td>Hepatitis B Virus DNA</td> <td>IU/mL</td> </tr> <tr> <td>85992623</td> <td>Hepatitis B Virus DNA</td> <td>Log IU/mL</td> </tr> </tbody> </table>		Result Code	Result Name	Unit of Measure	85997326	Hepatitis B Virus DNA	IU/mL	85992623	Hepatitis B Virus DNA	Log IU/mL
Result Code	Result Name	Unit of Measure									
85997326	Hepatitis B Virus DNA	IU/mL									
85992623	Hepatitis B Virus DNA	Log IU/mL									
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>2479</td> <td>Hepatitis B Virus MonitR, Chronic</td> </tr> <tr> <td>8137</td> <td>Hepatitis B Virus DNA, Quantitative, Real Time PCR</td> </tr> </tbody> </table>		Test Codes:	Name:	2479	Hepatitis B Virus MonitR, Chronic	8137	Hepatitis B Virus DNA, Quantitative, Real Time PCR			
Test Codes:	Name:										
2479	Hepatitis B Virus MonitR, Chronic										
8137	Hepatitis B Virus DNA, Quantitative, Real Time PCR										

AST (Aspartate Aminotransferase)															
Clinical Significance	AST is widely distributed throughout the tissues with significant amounts being in the heart and liver. Lesser amounts are found in skeletal muscles, kidneys, pancreas, spleen, lungs and brain. Injury to these tissues results in the release of the AST enzyme to general circulation. In myocardial infarction, serum AST may begin to rise within 6-8 hours after onset, peak within two days and return to normal by the fourth or fifth day post infarction. An increase in serum AST is also found with hepatitis, liver necrosis, cirrhosis and liver metastasis.														
Effective Date	3/10/2014														
Former Test Code	1345														
Test Code	822														
Instructions	Centrifuge serum or plasma specimens within 1 hour of collection, transfer serum or plasma to a sterile, plastic, screw-capped vial(s), and ship at room temperature.														
Reference Range	<table border="1"> <thead> <tr> <th>Age</th> <th>Male Reference Range</th> </tr> </thead> <tbody> <tr> <td>< 1 month:</td> <td>3- 51 U/L</td> </tr> <tr> <td>1-11 months:</td> <td>3- 65 U/L</td> </tr> <tr> <td>1-3 years:</td> <td>3- 56 U/L</td> </tr> <tr> <td>4-6 years:</td> <td>20-39 U/L</td> </tr> <tr> <td>7-19 years:</td> <td>12-32 U/L</td> </tr> <tr> <td>20-49 years:</td> <td>10-40 U/L</td> </tr> </tbody> </table>	Age	Male Reference Range	< 1 month:	3- 51 U/L	1-11 months:	3- 65 U/L	1-3 years:	3- 56 U/L	4-6 years:	20-39 U/L	7-19 years:	12-32 U/L	20-49 years:	10-40 U/L
Age	Male Reference Range														
< 1 month:	3- 51 U/L														
1-11 months:	3- 65 U/L														
1-3 years:	3- 56 U/L														
4-6 years:	20-39 U/L														
7-19 years:	12-32 U/L														
20-49 years:	10-40 U/L														

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	<table border="1"> <tr> <td>> or = 50 years:</td> <td>10-35 U/L</td> </tr> <tr> <td>Age</td> <td>Female Reference Range</td> </tr> <tr> <td>< 1 month:</td> <td>3- 49 U/L</td> </tr> <tr> <td>1-11 months:</td> <td>3- 79 U/L</td> </tr> <tr> <td>1-3 years:</td> <td>3- 69 U/L</td> </tr> <tr> <td>4-6 years:</td> <td>20-39 U/L</td> </tr> <tr> <td>7-19 years:</td> <td>12-32 U/L</td> </tr> <tr> <td>20-44 years:</td> <td>10-30 U/L</td> </tr> <tr> <td>> or = 45 years:</td> <td>10-35 U/L</td> </tr> </table>	> or = 50 years:	10-35 U/L	Age	Female Reference Range	< 1 month:	3- 49 U/L	1-11 months:	3- 79 U/L	1-3 years:	3- 69 U/L	4-6 years:	20-39 U/L	7-19 years:	12-32 U/L	20-44 years:	10-30 U/L	> or = 45 years:	10-35 U/L
> or = 50 years:	10-35 U/L																		
Age	Female Reference Range																		
< 1 month:	3- 49 U/L																		
1-11 months:	3- 79 U/L																		
1-3 years:	3- 69 U/L																		
4-6 years:	20-39 U/L																		
7-19 years:	12-32 U/L																		
20-44 years:	10-30 U/L																		
> or = 45 years:	10-35 U/L																		
Performing Site	Quest Diagnostics Nichols Institute, Valencia																		
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> <th>Unit of Measure</th> </tr> </thead> <tbody> <tr> <td>25002300</td> <td>Aspartate Aminotransferase</td> <td>U/L</td> </tr> </tbody> </table>	Result Code	Result Name	Unit of Measure	25002300	Aspartate Aminotransferase	U/L												
Result Code	Result Name	Unit of Measure																	
25002300	Aspartate Aminotransferase	U/L																	
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>5317</td> <td>Comprehensive Metabolic Panel</td> </tr> <tr> <td>5318</td> <td>Hepatic Function Panel</td> </tr> </tbody> </table>	Test Codes:	Name:	5317	Comprehensive Metabolic Panel	5318	Hepatic Function Panel												
Test Codes:	Name:																		
5317	Comprehensive Metabolic Panel																		
5318	Hepatic Function Panel																		

Cardio IQ™ Lipoprotein Fractionation, Ion Mobility																								
Effective Date	3/10/2014																							
Test Code	91604																							
Specimen Requirements	1 mL (0.25 mL minimum) serum Plasma is no longer acceptable																							
Reference Range	<table border="1"> <thead> <tr> <th>Cardio IQ™ Lipoprotein Fractionation, Ion Mobility</th> <th>Male</th> <th>Female</th> </tr> </thead> <tbody> <tr> <td>LDL Particle Number</td> <td>1016-2185 nmol/L</td> <td>1016-2185 nmol/L</td> </tr> <tr> <td>LDL Small</td> <td>123-441 nmol/L</td> <td>115-386 nmol/L</td> </tr> <tr> <td>LDL Medium</td> <td>167-465 nmol/L</td> <td>121-397 nmol/L</td> </tr> <tr> <td>HDL Large</td> <td>4334-10815</td> <td>5038-17886</td> </tr> <tr> <td>LDL Pattern</td> <td>A Pattern</td> <td>A Pattern</td> </tr> <tr> <td>LDL Peak Size</td> <td>> or = 218.2 Angstrom</td> <td>> or = 218.2 Angstrom</td> </tr> </tbody> </table>	Cardio IQ™ Lipoprotein Fractionation, Ion Mobility	Male	Female	LDL Particle Number	1016-2185 nmol/L	1016-2185 nmol/L	LDL Small	123-441 nmol/L	115-386 nmol/L	LDL Medium	167-465 nmol/L	121-397 nmol/L	HDL Large	4334-10815	5038-17886	LDL Pattern	A Pattern	A Pattern	LDL Peak Size	> or = 218.2 Angstrom	> or = 218.2 Angstrom		
Cardio IQ™ Lipoprotein Fractionation, Ion Mobility	Male	Female																						
LDL Particle Number	1016-2185 nmol/L	1016-2185 nmol/L																						
LDL Small	123-441 nmol/L	115-386 nmol/L																						
LDL Medium	167-465 nmol/L	121-397 nmol/L																						
HDL Large	4334-10815	5038-17886																						
LDL Pattern	A Pattern	A Pattern																						
LDL Peak Size	> or = 218.2 Angstrom	> or = 218.2 Angstrom																						
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano																							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> <th>Unit of Measure</th> </tr> </thead> <tbody> <tr> <td>86002760</td> <td>LDL Particle Number</td> <td>nmol/L</td> </tr> <tr> <td>86009431</td> <td>LDL Small</td> <td>nmol/L</td> </tr> </tbody> </table>	Result Code	Result Name	Unit of Measure	86002760	LDL Particle Number	nmol/L	86009431	LDL Small	nmol/L														
Result Code	Result Name	Unit of Measure																						
86002760	LDL Particle Number	nmol/L																						
86009431	LDL Small	nmol/L																						

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	86009433	LDL Medium	nmol/L
	86006295	HDL Large	nmol/L
	86002762	LDL Pattern	Pattern
	86002761	LDL Peak Size	Angstrom
Additional Information	Update report format		

Neutrophil Function, Oxidative Burst							
Message	**This test code is now available for New York patient testing at Quest Diagnostics Nichols Institute, San Juan Capistrano**						
Effective Date	3/10/2014						
Test Code	37523						
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano						
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>85992637</td> <td>%Oxidation Pos Neutrophils</td> </tr> <tr> <td>86010650</td> <td>Specimen Age</td> </tr> </tbody> </table>	Result Code	Result Name	85992637	%Oxidation Pos Neutrophils	86010650	Specimen Age
Result Code	Result Name						
85992637	%Oxidation Pos Neutrophils						
86010650	Specimen Age						
Additional Information	Update report format.						

Casein (f78) IgG					
Effective Date	3/17/2014				
Former Test Name	<i>Allergen - Casein IgG</i>				
Former Test Code	<i>F78G</i>				
Test Code	38100				
Specimen Requirements	0.3 mL (0.15 mL minimum) serum				
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days				
Set-up/Analytic Time	Set up: Sun, Wed, Fri; Report available: 1-3 days				
Reference Range	<2.0 mcg/mL				
Always Message	This test(s) was performed using a kit that has not been cleared or approved by the FDA. The analytical performance characteristics of this test have been determined by Quest Diagnostics Nichols Institute. This test, and any food specific allergen IgG result, should not be used for the diagnosis of allergic or atopic disease states (except for sensitivity to milk in neonates and gluten sensitivity). The use of food specific allergen IgG results should be restricted to the assessment of response to therapeutic interventions.				
Methodology	Immunoassay				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>85995765</td> <td>Casein (f78) IgG</td> </tr> </tbody> </table>	Result Code	Result Name	85995765	Casein (f78) IgG
Result Code	Result Name				
85995765	Casein (f78) IgG				

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Tests Affected	Test Codes:	Name:
	14788	Food Specific IgG Allergy (Adult) Panel
	14791	Food Specific IgG Allergy (Pediatric) Panel
	29496	Food Specific IgG Allergy Panel 1
	37334	Food Specific IgG Allergy Panel 2

Egg White (f1) IgG					
Effective Date	3/17/2014				
Former Test Name	Allergen - Egg White IgG				
Former Test Code	F1G				
Test Code	34678				
Specimen Requirements	0.3 mL (0.15 mL minimum) serum				
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days				
Set-up/Analytic Time	Set up: Sun, Wed, Fri; Report available: 1-3 days				
Reference Range	<2.0 mcg/mL				
Always Message	This test(s) was performed using a kit that has not been cleared or approved by the FDA. The analytical performance characteristics of this test have been determined by Quest Diagnostics Nichols Institute. This test, and any food specific allergen IgG result, should not be used for the diagnosis of allergic or atopic disease states (except for sensitivity to milk in neonates and gluten sensitivity). The use of food specific allergen IgG results should be restricted to the assessment of response to therapeutic interventions.				
Methodology	Immunoassay				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>85995766</td> <td>Egg White (f1) IgG</td> </tr> </tbody> </table>	Result Code	Result Name	85995766	Egg White (f1) IgG
Result Code	Result Name				
85995766	Egg White (f1) IgG				
Tests Affected	Test Codes:	Name:			
	14788	Food Specific IgG Allergy (Adult) Panel			
	14791	Food Specific IgG Allergy (Pediatric) Panel			
	29496	Food Specific IgG Allergy Panel 1			
	37334	Food Specific IgG Allergy Panel 2			

Latex (k82) IgE	
Effective Date	3/17/2014
Former Test Name	Allergen-Latex (Brazilian Rubber Tree) IgE

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Former Test Code	K82							
Test Code	8927							
Specimen Stability	Room temperature: 14 days Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: Daily; Report available: 1-3 day							
Reference Range	<0.35 kU/L							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55109005</td> <td>Latex (k82) IgE</td> </tr> <tr> <td>55109010</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	55109005	Latex (k82) IgE	55109010	Class
Result Code	Result Name							
55109005	Latex (k82) IgE							
55109010	Class							

Mace (f266) IgE								
Effective Date	3/17/2014							
Former Test Name	Allergen - Mace IgE							
Former Test Code	RF266							
Test Code	3025							
Specimen Requirements	0.3 mL (0.15 mL minimum) serum							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: daily; Report available: 1-14 days							
Reference Range	<0.35 kU/L							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55172905</td> <td>Mace (f266) IgE</td> </tr> <tr> <td>55172910</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	55172905	Mace (f266) IgE	55172910	Class
Result Code	Result Name							
55172905	Mace (f266) IgE							
55172910	Class							

Persimmon (Kaki Fruit) (f301) IgE		
Effective Date	3/17/2014	
Former Test Name	Allergen - Persimmon IgE	
Former Test Code	RF301	
Test Code	3252	

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Specimen Requirements	0.3 mL (0.15 mL minimum) serum							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: daily; Report available: 1-14 days							
Reference Range	<0.35 kU/L							
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, Valencia, CA. It has not been cleared or approved by the U.S. Food and Drug Administration. The FDA has determined that such clearance or approval is not necessary. Performance characteristics refer to the analytical performance of the test.							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55177305</td> <td>Persimmon (f301) IgE</td> </tr> <tr> <td>55177310</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	55177305	Persimmon (f301) IgE	55177310	Class
Result Code	Result Name							
55177305	Persimmon (f301) IgE							
55177310	Class							

Salt Grass (g203) IgE								
Effective Date	3/17/2014							
Former Test Name	Allergen - Salt Grass IgE							
Former Test Code	RG203							
Test Code	2350							
Specimen Requirements	0.3 mL (0.15 mL minimum) serum							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: daily; Report available: 1-14 days							
Reference Range	<0.35 kU/L							
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, Valencia, CA. It has not been cleared or approved by the U.S. Food and Drug Administration. The FDA has determined that such clearance or approval is not necessary. Performance characteristics refer to the analytical performance of the test.							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55102205</td> <td>Salt Grass (g203) IgE</td> </tr> <tr> <td>55102210</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	55102205	Salt Grass (g203) IgE	55102210	Class
Result Code	Result Name							
55102205	Salt Grass (g203) IgE							
55102210	Class							

Staphylococcal Enterotoxin A (m80) IgE

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Effective Date	3/17/2014							
<i>Former Test Name</i>	<i>Allergen - Staphylococcal Enterotoxin A IgE</i>							
<i>Former Test Code</i>	<i>M80</i>							
Test Code	91996							
Specimen Requirements	0.3 mL (0.15 mL minimum) serum							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: daily; Report available: 1-14 days							
Reference Range	<0.35 kU/L							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86010458</td> <td>S. Enterotoxin A (m80) IgE</td> </tr> <tr> <td>86010485</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	86010458	S. Enterotoxin A (m80) IgE	86010485	Class
Result Code	Result Name							
86010458	S. Enterotoxin A (m80) IgE							
86010485	Class							

Staphylococcal Enterotoxin B (m81) IgE								
Effective Date	3/17/2014							
<i>Former Test Name</i>	<i>Allergen - Staphylococcal Enterotoxin B IgE</i>							
<i>Former Test Code</i>	<i>M81</i>							
Test Code	91997							
Specimen Requirements	0.3 mL (0.15 mL minimum) serum							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: daily; Report available: 1-14 days							
Reference Range	<0.35 kU/L							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86010459</td> <td>S. Enterotoxin B (m81) IgE</td> </tr> <tr> <td>86010633</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	86010459	S. Enterotoxin B (m81) IgE	86010633	Class
Result Code	Result Name							
86010459	S. Enterotoxin B (m81) IgE							
86010633	Class							

Ulocladium chartarum (m204) IgE		
Effective Date	3/17/2014	

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Former Test Name	Allergen - Ulocladium chartarum IgE							
Former Test Code	RM204							
Test Code	3412							
Specimen Requirements	0.3 mL (0.15 mL minimum) serum							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days							
Set-up/Analytic Time	Set up: daily; Report available: 1-14 days							
Reference Range	<0.35 kU/L							
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, Valencia, CA. It has not been cleared or approved by the U.S. Food and Drug Administration. The FDA has determined that such clearance or approval is not necessary. Performance characteristics refer to the analytical performance of the test.							
Methodology	Immunoassay							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55156405</td> <td>Ulocladium chartarum (m204) IgE</td> </tr> <tr> <td>55156410</td> <td>Class</td> </tr> </tbody> </table>		Result Code	Result Name	55156405	Ulocladium chartarum (m204) IgE	55156410	Class
Result Code	Result Name							
55156405	Ulocladium chartarum (m204) IgE							
55156410	Class							

Wheat (f4) IgG						
Effective Date	3/17/2014					
Former Test Name	Allergen - Wheat IgG					
Former Test Code	F4G					
Test Code	30337					
Specimen Requirements	0.3 mL (0.15 mL minimum) serum					
Specimen Stability	Room temperature and Refrigerated: 14 days Frozen: 30 days					
Set-up/Analytic Time	Set up: Sun, Wed, Fri; Report available: 1-3 days					
Reference Range	<2.0 mcg/mL					
Always Message	This test(s) was performed using a kit that has not been cleared or approved by the FDA. The analytical performance characteristics of this test have been determined by Quest Diagnostics Nichols Institute. This test, and any food specific allergen IgG result, should not be used for the diagnosis of allergic or atopic disease states (except for sensitivity to milk in neonates and gluten sensitivity). The use of food specific allergen IgG results should be restricted to the assessment of response to therapeutic interventions.					
Methodology	Immunoassay					
Performing Site	Quest Diagnostics Nichols Institute, Valencia					
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>85995770</td> <td>Wheat (f4) IgG</td> </tr> </tbody> </table>		Result Code	Result Name	85995770	Wheat (f4) IgG
Result Code	Result Name					
85995770	Wheat (f4) IgG					

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Tests Affected	Test Codes:	Name:
	14788	Food Specific IgG Allergy (Adult) Panel
	14791	Food Specific IgG Allergy (Pediatric) Panel
	29496	Food Specific IgG Allergy Panel 1
	37334	Food Specific IgG Allergy Panel 2

QuantiFERON®-TB Gold, (Draw Site Incubated)	
Clinical Significance	QuantiFERON® TB Gold IT is an indirect test for <i>M. tuberculosis</i> infection (including disease) and is intended for use in conjunction with risk assessment, radiography and other medical and diagnostic evaluations.
Effective Date	3/31/2014
Former Test Code	6110
Test Code	16603
CPT Codes	86480
Specimen Requirements	1 mL (0.8 mL minimum) whole blood collected in each QTF-Nil control, QTF-TB AG and QTF-Mitogen control tube
Reject Criteria	Received >72 hours (room temperature); received frozen; specimens not collected in QuantiFERON® Gold IT blood collection tubes Remove: >72 hours (pre-incubated)
Instructions	<p>1. For each patient, collect 1 mL of blood by venipuncture directly into each of the three (3) unique QuantiFERON®- TB Gold IT blood collection tubes. Tubes must be at room temperature prior to collection. Under or overfilling of the tubes may lead to erroneous results.</p> <p>2. Shake them ten (10) times just firmly enough to ensure the entire surface of the tube is coated with blood, to solubilize antigens on tube walls.</p> <p>3. Incubate the three (3) tubes upright at 36-38° C within 16 hours for 16 to 24 hours. The specimen rack should sit off the floor of the incubator. If the blood is not incubated immediately after collection, remixing of the tubes by inverting 10 times must be repeated immediately prior to incubation.</p> <p>4. Following incubation either:</p> <p>A. Immediately transport the three (3) transport tubes to Quest Diagnostics between 4° and 27° C. Samples will be stable for 72 hours at 4°-27° C (room temperature or refrigerated).</p> <p>OR</p> <p>B. Centrifuge each of the three (3) incubated collection tubes for 15 minutes at 2000 to 3000 RCF (g). Label with patient name, identification number, and date of collection. Deliver to Quest Diagnostics at 2-8 degrees C. Samples will be stable for 28 days at 2-8 degrees C (refrigerated).</p>
Transport Temperature	See instructions
Specimen Stability	<p>Uncentrifuged specimens: Room temperature and Refrigerated: 72 hours Frozen: Unacceptable</p> <p>Centrifuged specimens: Room temperature: 72 hours Refrigerated: 28 days Frozen: Unacceptable</p>
Set-up/Analytic Time	Set up: Daily; Report available: 1-2 days
Reference Range	Negative
Always Message	The Nil tube value is used to determine if the patient has a preexisting immune response which could cause a false-positive reading on the test. In order for a test to be valid, the Nil tube must have a value of

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	<p>less than or equal to 8.0 IU/mL.</p> <p>The mitogen control tube is used to assure the patient has a healthy immune status and also serves as a control for correct blood handling and incubation. It is used to detect false-negative readings. The mitogen tube must have a gamma interferon value of greater than or equal to 0.5 IU/mL higher than the value of the Nil tube.</p> <p>The TB antigen tube is coated with the M. tuberculosis specific antigens. For a test to be considered positive, the TB antigen tube value minus the Nil tube value must be greater than or equal to 0.35 IU/mL.</p> <p>For additional information, please refer to http://education.questdiagnostics.com/faq/QFT (This link is being provided for informational/educational purposes only.)</p> <p>Data on the performance of the test in children younger than 5 years of age are limited, and the CDC advises that caution is warranted when using the assay in children aged <5 years (MMWR 2010;59(RR-05):1-25).</p>															
Methodology	Immunoassay															
Performing Site	Quest Diagnostics Nichols Institute, Valencia															
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> <th>Unit of Measure</th> </tr> </thead> <tbody> <tr> <td>40023705</td> <td>QuantiFERON(R)TB,Incubated</td> <td></td> </tr> <tr> <td>40023710</td> <td>NIL</td> <td>IU/mL</td> </tr> <tr> <td>40023730</td> <td>Mitogen-NIL</td> <td>IU/mL</td> </tr> <tr> <td>40023720</td> <td>TB-NIL</td> <td>IU/mL</td> </tr> </tbody> </table>	Result Code	Result Name	Unit of Measure	40023705	QuantiFERON(R)TB,Incubated		40023710	NIL	IU/mL	40023730	Mitogen-NIL	IU/mL	40023720	TB-NIL	IU/mL
Result Code	Result Name	Unit of Measure														
40023705	QuantiFERON(R)TB,Incubated															
40023710	NIL	IU/mL														
40023730	Mitogen-NIL	IU/mL														
40023720	TB-NIL	IU/mL														

Discontinued Tests

<i>Borrelia burgdorferi</i> IgG & IgM Abs [EIA] Plus C6 Peptide	
Effective Date	3/10/2014
Test Code	8954
Additional Information	The recommended alternative is test code 8941- <i>Borrelia burgdorferi</i> IgG & IgM Abs [EIA]

<i>Borrelia burgdorferi</i> IgG & IgM Abs + C6 Peptide w/Rfx IB [CDC]	
Effective Date	3/10/2014
Test Code	8956
Additional Information	<p>The recommended alternatives are test codes:</p> <ul style="list-style-type: none"> • 8942-<i>Borrelia burgdorferi</i> IgG & IgM Abs w/Reflex IB + Bands [CDC] • 8942NY-<i>Borrelia burgdorferi</i> IgG & IgM Abs w/Reflex IB + Bands [NY]

Heparin Anti-Xa (Unfractionated Heparin)	
Effective Date	3/10/2014
Test Code	S52557
Additional Information	The recommended alternative is 30292 - Heparin, Anti-Xa, in the Test Changes Section.

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Lyme Disease C6 Antibodies, Total, ELISA					
Effective Date	3/10/2014				
Test Code	8944				
Additional Information	There is no recommended alternative.				
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>8946</td> <td>Borrelia Burgdorferi C6 Peptide Ab DetectR w/Rfx IB [CDC]</td> </tr> </tbody> </table>	Test Codes:	Name:	8946	Borrelia Burgdorferi C6 Peptide Ab DetectR w/Rfx IB [CDC]
Test Codes:	Name:				
8946	Borrelia Burgdorferi C6 Peptide Ab DetectR w/Rfx IB [CDC]				

Magnesium RBC [623Z][NY]	
Effective Date	3/10/2014
Test Code	S49785NY
Additional Information	The recommended alternative is test code 623 - Magnesium, RBC

Allergen - Latex Enhanced IgE	
Effective Date	3/17/2014
Test Code	RK82
Additional Information	The recommended alternative is 8927 -Latex (f82) IgE.

Test Send Outs (Referrals)

Due to these assays being performed by outside vendors, we are unable to use our normal method of communication. Some of the changes listed in this document may be effective in less than 30 days. Please note the individual effective dates below, as these changes may require **IMMEDIATE ACTION**.

Amphetamines Panel, Serum/Plasma	
Effective Date	3/3/2014
Test Code	S47280
Specimen Requirements	<p>Preferred: 1 mL (0.4 mL minimum) serum collected in a red-top tube (no gel)</p> <p>Acceptable: Plasma collected in an EDTA (lavender-top) or EDTA (pink-top) tube</p>

Fluphenazine, Serum/Plasma	
Effective Date	3/3/2014
Test Code	92021
Specimen Requirements	<p>Preferred: 2 mL (0.7 mL minimum) serum collected in a red-top tube (no gel)</p> <p>Acceptable: Plasma collected in an EDTA (lavender-top or pink-top) tube</p> <p>Serum collected in a no additive (blue-top) tube or plasma collected in an EDTA (royal-blue-top) tube is not acceptable.</p>

February 2014 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Reject Criteria	Received room temperature; polymer gel separation tube					
Transport Temperature	Refrigerated					
Specimen Stability	Room temperature: 24 hours Refrigerated and Frozen: 14 days					
Set-up/Analytic Time	Set up: Wed; Report available: 3 days					
Reference Range	Schizophrenic patients maintained with depot injections of fluphenazine decanoate had the following plasma fluphenazine concentrations: 1 to 3 ng/mL following 12.5 mg per week 4 to 7 ng/mL following 25 mg per week 5 to 17 ng/mL following 50 mg per week Healthy subjects given single oral doses of 5 mg fluphenazine had peak plasma concentrations averaging 0.6 ng/mL (SEM +/- 0.1 ng/mL)					
Units Of Measure	ng/mL					
Methodology	High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)					
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86010491</td> <td>Fluphenazine, Serum/Plasma</td> </tr> </tbody> </table>		Result Code	Result Name	86010491	Fluphenazine, Serum/Plasma
Result Code	Result Name					
86010491	Fluphenazine, Serum/Plasma					

Fluphenazine, Serum/Plasma	
Effective Date	3/3/2014
Test Code	S41755
Additional Information	This test will be discontinued and the recommended alternative is test code 92021- Fluphenazine, Serum/Plasma.