

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

ANNOUNCEMENTS		
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.		
Test Code	Test Name	Page #
	HLA Typing	2

NEW TESTS			
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.			
Test Code	Test Name	Effective Date	Page #
91427	Melanoma, Chromosomal Microarray, ClariSure® Oligo-SNP	4/9/2013	3
91136	Secretory Phospholipase A2, Type IIA (sPLA2)	4/9/2013	4
91461	Lynch Syndrome Panel	4/16/2013	4
91584	Lynch Syndrome, MLH1 Familial Deletion/Duplication	4/16/2013	5
91460	Lynch Syndrome, MLH1 Sequencing and Deletion/Duplication	4/16/2013	6
91459	Lynch Syndrome, MSH2 Familial Deletion/Duplication (Including EPCAM)	4/16/2013	6
91471	Lynch Syndrome, MSH2 Sequencing and Deletion/Duplication (Including EPCAM)	4/16/2013	7
91458	Lynch Syndrome, MSH6 Sequencing and Deletion/Duplication	4/16/2013	8
91463	Lynch Syndrome, PMS2 Familial Deletion/Duplication	4/16/2013	9
91457	Lynch Syndrome, PMS2 Sequencing and Deletion/Duplication	4/16/2013	9
17303	Chlamydia trachomatis DNA, SDA	4/30/2013	10
17615	Chlamydia trachomatis DNA, SDA, Pap Vial	4/30/2013	11
17618	Chlamydia/N. gonorrhoeae DNA, SDA, Pap Vial	4/30/2013	12
17305	Chlamydia/N.gonorrhoeae DNA, SDA	4/30/2013	13
17304	Neisseria gonorrhoeae DNA, SDA	4/30/2013	14
17617	Neisseria gonorrhoeae DNA, SDA, Pap Vial	4/30/2013	15

TEST CHANGES				
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.				
Test Code	Former Test Code	Test Name	Effective Date	Page #
90521		Trichomonas vaginalis RNA, Qualitative TMA, Pap Vial	3/11/2013	15
10051		Hepatitis C Viral RNA, Quantitative Real-Time PCR w/Rfx to Qualitative TMA	3/26/2013	16
8137		Hepatitis B Virus DNA, Quantitative, Real-Time PCR	4/1/2013	16
7577		Hepatitis C Viral RNA, Quantitative, Real-Time PCR	4/1/2013	17
7485A		HIV-1 RNA, Quantitative, Real-Time PCR	4/1/2013	18
36712		Sirolimus, LC/MS/MS	4/1/2013	18
70007		Tacrolimus, Highly Sensitive, LC/MS/MS	4/1/2013	19

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

S51345		21-Hydroxylase Antibody [37916X]	4/8/2013	19
6399		CBC (Includes Differential and Platelets)	4/8/2013	19
S51491		Dihydrotestosterone, LC/MS/MS [90567]	4/8/2013	21
S49575		Estriol, LC/MS/MS, Serum [34883]	4/8/2013	21
S51290		Muramidase (Lysozyme), Serum	4/8/2013	21
S51657		Porphobilinogen, Quantitative, 24-Hour Urine	4/8/2013	21
S51658		Porphobilinogen, Quantitative, Random Urine	4/8/2013	22
S48558		Porphyrins, Fractionated, Quantitative, Random Urine	4/8/2013	22
S51406		Porphyrins, Total, Plasma	4/8/2013	22
4948		Protein C Antigen	4/8/2013	23
3060		Thyroid Abs Evaluation	4/16/2013	23
8761		Varicella-Zoster Virus Antibody (IgG)	4/30/2013	23

REDIRECTS

Please Note: Not all test codes assigned to each assay are listed in the table of contents.
Please refer to the complete listing on the page numbers indicated.

Test Code	Former Test Code	Test Name	Effective Date	Page #
30551	1092	TSI (Thyroid Stimulating Immunoglobulin)	4/16/2013	24
5042X	3515W	Vitamin B1 (Thiamine), Blood, LC/MS/MS	4/23/2013	24
16500	2040	HE4, Ovarian Cancer Monitoring	5/6/2013	25

DISCONTINUED TESTS

Please Note: Not all test codes assigned to each assay are listed in the table of contents.
Please refer to the complete listing on the page numbers indicated.

Test Code	Test Name	Effective Date	Page #
2445	Hepatitis C Antibody w/Reflex RIBA	4/16/2013	26
1091	Thyroid Stimulating Immunoglobulins w/TSH	4/16/2013	26
S51330	Chlamydia/N. gonorrhoeae DNA, SDA	4/30/2013	26

NY UPDATE

Please Note: Not all test codes assigned to each assay are listed in the table of contents.
Please refer to the complete listing on the page numbers indicated.

Test Code	Test Name	Page #
90569	Herpes Simplex Virus, Type 1/2 DNA, Real-Time PCR, PAP Vial	26

Announcements

HLA Typing	
Effective Date	3/11/2013

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Reference Range	<p>The HLA DNA Typing results will change to include a colon (:) as mandated by the WHO Nomenclature Committee.</p> <p>Example: for Test Code 17397X [17397] - HLA-A High Resolution SBT Typing: HLA-A* 02:01 HLA-A* 03:01</p> <p>Additional details may be found in http://hla.alleles.org/nomenclature/nomenc_reports.html : <i>Nomenclature for factors of the HLA system, 2010</i>, Tissue Antigens 75, 291-455.</p>											
Performing Site	Quest Diagnostics Nichols Institute, Chantilly											
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>S51547</td> <td>HLA-A High Resolution SBT Typing [17397X]</td> </tr> <tr> <td>S51548</td> <td>HLA-B High Resolution SBT Typing [17396X]</td> </tr> <tr> <td>S52488</td> <td>HLA-DRB1 High Resolution SBT Typing [17393X]</td> </tr> <tr> <td>S51409</td> <td>HLA-DQB1 High Resolution SBT Typing [17394X]</td> </tr> </tbody> </table>		Test Codes:	Name:	S51547	HLA-A High Resolution SBT Typing [17397X]	S51548	HLA-B High Resolution SBT Typing [17396X]	S52488	HLA-DRB1 High Resolution SBT Typing [17393X]	S51409	HLA-DQB1 High Resolution SBT Typing [17394X]
Test Codes:	Name:											
S51547	HLA-A High Resolution SBT Typing [17397X]											
S51548	HLA-B High Resolution SBT Typing [17396X]											
S52488	HLA-DRB1 High Resolution SBT Typing [17393X]											
S51409	HLA-DQB1 High Resolution SBT Typing [17394X]											

New Test Offerings

The following tests will be available through Quest Diagnostics on the dates indicated below.

Melanoma, Chromosomal Microarray, ClariSure® Oligo-SNP									
Clinical Significance	This assay detects alterations in copy number or allelic imbalance. It can be used in the differential diagnosis of melanocytic lesions (e.g., classifying indeterminate lesions as benign or malignant). The Illumina Infinium oligo-SNP array reveals copy number variations and loss of heterozygosity (LOH) across the whole genome. The allelic copy number alterations, often complex, are associated with malignant melanoma. The assay is performed using DNA extracted from FFPE samples.								
Effective Date	4/9/2013								
Test Code	91427								
CPT Codes	81406								
Specimen Requirements	Formalin fixed paraffin embedded tissue block								
Instructions	<p>Formalin fixed paraffin-embedded block. Specimen must be fixed in 10% neutral buffered formalin. Fixation between 6 and 48 hours is recommended. Pathology report must accompany paraffin block. Information required in this report include: Physician identification, specimen identifiers (case and block number), specimen site and type, tissue processing used (routine or microwave), type of fixative, time and duration of fixation, and pathologic diagnosis.</p> <p>Ship at room temperature. Do not freeze.</p>								
Transport Temperature	Room temperature								
Specimen Stability	See Instruction								
Set-up/Analytic Time	Set up: As needed; Report available: within 22 days								
Methodology	Microarray								
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano								
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Type</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009049</td> <td></td> <td>Melanoma,ClariSure, Oligo</td> </tr> </tbody> </table>			Result Code	Type	Result Name	86009049		Melanoma,ClariSure, Oligo
Result Code	Type	Result Name							
86009049		Melanoma,ClariSure, Oligo							

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	86001293	AOE	Specimen Source:
	86007540	AOE	Diagnosis:
	86007541	AOE	Physician Name:
	86007542	AOE	Physician Phone #:
	86007543	AOE	Physician NPI

Secretary Phospholipase A2, Type IIA (sPLA2)					
Clinical Significance	Higher levels of sPLA2 (individually or combined with other risk factors) can be used for risk stratification in prediction of atherosclerosis/CVD and estimation for recurrent events after ACS.				
Effective Date	4/9/2013				
Test Code	91136				
CPT Codes	83520				
Specimen Requirements	2 mL plasma collected in an EDTA (lavender-top) tube				
Reject Criteria	Serum; heparin; other plasma samples; grossly icteric				
Instructions	Separate plasma from cells as soon as possible after centrifugation. Transfer plasma to a Polypropylene transport tube for shipping.				
Transport Temperature	Room temperature				
Specimen Stability	Room temperature and Refrigerated: 7 days Frozen: 21 days Frozen -70 Degrees: 28 days				
Set-up/Analytic Time	Set up: Wed; Report available: 7 days				
Reference Range	Males: < or = 4.56 ng/mL Females: < or = 6.28 ng/mL				
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.				
Methodology	Enzyme Immunoassay				
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86008762</td> <td>sPLA2</td> </tr> </tbody> </table>	Result Code	Result Name	86008762	sPLA2
Result Code	Result Name				
86008762	sPLA2				

Lynch Syndrome Panel	
Clinical Significance	Germline mutations in DNA mismatch repair genes are the major cause of Hereditary Nonpolyposis Colorectal Cancer (HNPCC), also called Lynch Syndrome. HNPCC is typically inherited as autosomal dominant. The major genes involved in HNPCC are MLH1, MSH2, MSH6, and PMS2. HNPCC accounts for between 5% and 10% of colon cancers. Sequence analysis and deletion/duplication analysis of MLH1, MSH2, MSH6, and PMS2 are performed in this panel.
Effective Date	4/16/2013
Test Code	91461
CPT Codes	81295, 81297, 81292, 81294, 81298, 81300, 81317, 81319, 81403

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Specimen Requirements	4 mL (2 mL minimum) whole blood collected in an EDTA (lavender-top) tube [x5]																					
Reject Criteria	Samples received frozen																					
Instructions	Whole blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.																					
Transport Temperature	Room temperature																					
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable																					
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days																					
Reference Range	Accompanies report																					
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.																					
Methodology	Sequencing, Multiplex Ligation-dependent Probe Amplification (MLPA), Semi-quantitative Fluorescent Polymerase Chain Reaction																					
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano																					
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009141</td> <td>Summary</td> </tr> <tr> <td>86009308</td> <td>MSH2 Sequencing</td> </tr> <tr> <td>86009137</td> <td>MSH2/EPCAM Del/Dup</td> </tr> <tr> <td>86009307</td> <td>MLH1 Sequencing</td> </tr> <tr> <td>86009138</td> <td>MLH1 Deletion, Duplication</td> </tr> <tr> <td>86009306</td> <td>MSH6 Sequencing</td> </tr> <tr> <td>86008784</td> <td>MSH6 Deletion, Duplication</td> </tr> <tr> <td>86009135</td> <td>PMS2 Sequencing</td> </tr> <tr> <td>86009136</td> <td>PMS2 Deletion, Duplication</td> </tr> </tbody> </table>		Result Code	Result Name	86009141	Summary	86009308	MSH2 Sequencing	86009137	MSH2/EPCAM Del/Dup	86009307	MLH1 Sequencing	86009138	MLH1 Deletion, Duplication	86009306	MSH6 Sequencing	86008784	MSH6 Deletion, Duplication	86009135	PMS2 Sequencing	86009136	PMS2 Deletion, Duplication
Result Code	Result Name																					
86009141	Summary																					
86009308	MSH2 Sequencing																					
86009137	MSH2/EPCAM Del/Dup																					
86009307	MLH1 Sequencing																					
86009138	MLH1 Deletion, Duplication																					
86009306	MSH6 Sequencing																					
86008784	MSH6 Deletion, Duplication																					
86009135	PMS2 Sequencing																					
86009136	PMS2 Deletion, Duplication																					

Lynch Syndrome, MLH1 Familial Deletion/Duplication	
Clinical Significance	This assay detects deletions and/or duplication involving one or more exons, as well as the entire MLH1 gene. This assay cannot detect smaller genetic alterations, such as point mutations affecting amino acid coding or mRNA splicing. The assay cannot detect mutations affecting MLH1 gene regions not examined in the assay (including most of the intronic regions). This assay cannot rule out the presence of germline mutations involving other mismatch repair genes.
Effective Date	4/16/2013
Test Code	91584
CPT Codes	81294
Specimen Requirements	4 mL whole blood collected in an EDTA (lavender-top) tube
Reject Criteria	Samples received frozen
Instructions	Whole blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.
Transport Temperature	Room temperature

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable				
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days				
Reference Range	Accompanies report				
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.				
Methodology	Semi-quantitative Fluorescent Polymerase Chain Reaction				
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009138</td> <td>MLH1 Deletion, Duplication</td> </tr> </tbody> </table>	Result Code	Result Name	86009138	MLH1 Deletion, Duplication
Result Code	Result Name				
86009138	MLH1 Deletion, Duplication				

Lynch Syndrome, MLH1 Sequencing and Deletion/Duplication			
Clinical Significance	<p>Sequencing: This test should be offered to patients with colorectal cancer who meet the Bethesda criteria established by the National Cancer Institute (NCI), and/or that their tumors display microsatellite instability (MSI) or replication error (RER) phenotype. Detection of a germline mutation in one or the mismatch repair genes helps to establish a clinical diagnosis of HNPCC in affected patients. Mutation detection in at risk family member allows predictive diagnosis of the disease and thus intensive screening and early intervention of cancer.</p> <p>Deletion/Duplication: This assay detects deletions and/or duplication involving one or more exons, as well as the entire MLH1 gene. This assay cannot detect smaller genetic alterations, such as point mutations affecting amino acid coding or mRNA splicing. The assay cannot detect mutations affecting MLH1 gene regions not examined in the assay (including most of the intronic regions). This assay cannot rule out the presence of germline mutations involving other mismatch repair genes.</p>		
Effective Date	4/16/2013		
Test Code	91460		
CPT Codes	81292, 81294		
Specimen Requirements	4 mL (2 mL minimum) whole blood collected in an EDTA (lavender-top) tube [x2]		
Reject Criteria	Samples received frozen		
Instructions	Whole Blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.		
Transport Temperature	Room temperature		
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable		
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days		
Reference Range	Accompanies report		
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.		
Methodology	Polymerase Chain Reaction and Dye-Terminator Sequencing Reaction; Semi-quantitative Fluorescent Polymerase Chain Reaction		
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano		
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> </tbody> </table>	Result Code	Result Name
Result Code	Result Name		

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	86009307	MLH1 Sequencing
	86009138	MLH1 Deletion, Duplication

Lynch Syndrome, MSH2 Familial Deletion/Duplication (Including EPCAM)					
Clinical Significance	This assay detects deletions and/or duplications involving one of more exons for the MSH2 gene, and exons 8 and 9 for the EPCAM gene. This assay cannot detect smaller genetic alterations, such as point mutations affecting amino acid coding or mRNA splicing. This assay cannot detect mutations affecting MSH2 and EPCAM gene regions not examined in the assay (including most of the intronic regions). This assay cannot rule out the presence of germline mutations involving other mismatch repair genes.				
Effective Date	4/16/2013				
Test Code	91459				
CPT Codes	81297, 81403				
Specimen Requirements	5 mL whole blood collected in an EDTA (lavender-top) tube				
Reject Criteria	Samples received frozen				
Instructions	Whole Blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do Not Freeze.				
Transport Temperature	Room temperature				
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable				
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days				
Reference Range	Accompanies report				
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.				
Methodology	Semi-quantitative Fluorescent Polymerase Chain Reaction				
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009137</td> <td>MSH2/EPCAM Del/Dup</td> </tr> </tbody> </table>	Result Code	Result Name	86009137	MSH2/EPCAM Del/Dup
Result Code	Result Name				
86009137	MSH2/EPCAM Del/Dup				

Lynch Syndrome, MSH2 Sequencing and Deletion/Duplication (Including EPCAM)	
Clinical Significance	<p>Sequencing: This test should be offered to patients with colorectal cancer who meet the Bethesda criteria established by the National Cancer Institute (NCI), and/or that their tumors display microsatellite instability (MSI) or replication error (RER) phenotype. Detection of a germline mutation in one or the mismatch repair genes helps to establish a clinical diagnosis of HNPCC in affected patients. Mutation detection in at risk family member allows predictive diagnosis of the disease and thus intensive screening and early intervention of cancer.</p> <p>Deletion/Duplication: This assay detects deletions and/or duplication involving one or more exons for the MSH2 gene, and exons 8 and 9 for the EPCAM gene. This assay cannot detect smaller genetic alterations, such as point mutations affecting amino acid coding or mRNA splicing. This assay cannot detect mutations affecting MSH2 and EPCAM gene regions not examined in the assay (including most of the intronic regions). This assay cannot rule out the presence of germline mutations involving other mismatch repair genes.</p>
Effective Date	4/16/2013
Test Code	91471

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

CPT Codes	81295, 81297, 81403							
Specimen Requirements	4 mL (2 mL minimum) whole blood collected in an EDTA (lavender-top) tube [x2]							
Reject Criteria	Samples received frozen							
Instructions	Whole Blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable							
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days							
Reference Range	Accompanies report							
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.							
Methodology	Polymerase Chain Reaction and Dye-Terminator Sequencing Reaction, Semi-quantitative Fluorescent Polymerase Chain Reaction							
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009308</td> <td>MSH2 Sequencing</td> </tr> <tr> <td>86009137</td> <td>MSH2/EPCAM Del/Dup</td> </tr> </tbody> </table>		Result Code	Result Name	86009308	MSH2 Sequencing	86009137	MSH2/EPCAM Del/Dup
Result Code	Result Name							
86009308	MSH2 Sequencing							
86009137	MSH2/EPCAM Del/Dup							

Lynch Syndrome, MSH6 Sequencing and Deletion/Duplication	
Clinical Significance	<p>Sequencing: This test should be offered to patients with colorectal cancer who meet the Bethesda criteria established by the National Cancer Institute (NCI), and/or that their tumors display microsatellite instability (MSI) or replication error (RER) phenotype. Detection of a germline mutation in one or the mismatch repair genes helps to establish a clinical diagnosis of HNPCC in affected patients. Mutation detection in at risk family member allows predictive diagnosis of the disease and thus intensive screening and early intervention of cancer.</p> <p>Deletion/Duplication: This assay detects deletion and/or duplication involving one or more exons, as well as the entire MSH6 gene. This assay cannot detect smaller alterations, such as point mutations affecting amino acid coding or mRNA splicing. The assay cannot detect mutations affecting gene regions not examined in the assay (including most of the intronic regions). This assay cannot rule out the presence of germline mutations involving other mismatch repair genes.</p>
Effective Date	4/16/2013
Test Code	91458
CPT Codes	81298, 81300
Specimen Requirements	4 mL (2mL minimum) whole blood collected in an EDTA (lavender-top) tube [x2]
Reject Criteria	Samples received frozen
Instructions	Whole Blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.
Transport Temperature	Room temperature
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Reference Range	Accompanies report							
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.							
Methodology	Polymerase Chain Reaction, Dye Terminator Sequencing Reaction and Semi-quantitative Fluorescent Polymerase Chain Reaction							
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009306</td> <td>MSH6 Sequencing</td> </tr> <tr> <td>86008784</td> <td>MSH6 Deletion, Duplication</td> </tr> </tbody> </table>		Result Code	Result Name	86009306	MSH6 Sequencing	86008784	MSH6 Deletion, Duplication
Result Code	Result Name							
86009306	MSH6 Sequencing							
86008784	MSH6 Deletion, Duplication							

Lynch Syndrome, PMS2 Familial Deletion/Duplication					
Clinical Significance	Hereditary non-polyposis colorectal cancer (HNPCC) is an autosomal dominant genetic disorder characterized by the familial accumulation of early onset colorectal, endometrial and other tumors, and accounts for approximately 3-5% of all colon cancer cases. HNPCC is caused by germline mutations in four genes involved in the DNA mismatch repair pathway: MLH1, MSH2, MSH6, and PMS2. MLH1 and MSH2 mutations account for approximately 90% of HNPCC cases, MSH6 accounts for approximately 7-10%, while PMS2 accounts for less than 5%. The prevalence of PMS2 mutations is not well established but some studies have shown that large deletions account for approximately 21-37% of PMS2 mutations.				
Effective Date	4/16/2013				
Test Code	91463				
CPT Codes	81319				
Specimen Requirements	4 mL whole blood collected in an EDTA (lavender-top) tube				
Reject Criteria	Samples received frozen				
Instructions	Whole blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.				
Transport Temperature	Room temperature				
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable				
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days				
Reference Range	Accompanies report				
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.				
Methodology	Multiplex Ligation-dependent Probe Amplification (MLPA)				
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009136</td> <td>PMS2 Deletion, Duplication</td> </tr> </tbody> </table>	Result Code	Result Name	86009136	PMS2 Deletion, Duplication
Result Code	Result Name				
86009136	PMS2 Deletion, Duplication				

Lynch Syndrome, PMS2 Sequencing and Deletion/Duplication

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Clinical Significance	Hereditary non-polyposis colorectal cancer (HNPCC) is an autosomal dominant genetic disorder characterized by the familial accumulation of early onset colorectal, endometrial and other tumors, and accounts for approximately 3-5% of all colon cancer cases. HNPCC is caused by germline mutations in four genes involved in the DNA mismatch repair pathway: MLH1, MSH2, MSH6, and PMS2. MLH1 and MSH2 mutations account for approximately 90% of HNPCC cases, MSH6 accounts for approximately 7-10%, while PMS2 accounts for less than 5%. The prevalence of PMS2 mutations is not well established but some studies have shown that large deletions account for approximately 21-37% of PMS2 mutations. The prevalence of sequence variants is not well established.							
Effective Date	4/16/2013							
Test Code	91457							
CPT Codes	81317, 81319							
Specimen Requirements	4 mL whole blood collected in an EDTA (lavender-top) tube							
Reject Criteria	Samples received frozen							
Instructions	Whole blood: Normal phlebotomy procedure. Specimen stability is crucial. Store and ship ambient immediately. Do not freeze.							
Transport Temperature	Room temperature							
Specimen Stability	Room temperature and Refrigerated: 8 days Frozen: Unacceptable							
Set-up/Analytic Time	Set up: Thurs; Report available: 15 days							
Reference Range	Accompanies report							
Always Message	This test was developed and its performance characteristics have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. Performance characteristics refer to the analytical performance of the test.							
Methodology	DNA Sequencing and Multiplex Ligation-dependent Probe Amplification (MLPA)							
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86009135</td> <td>PMS2 Sequencing</td> </tr> <tr> <td>86009136</td> <td>PMS2 Deletion, Duplication</td> </tr> </tbody> </table>		Result Code	Result Name	86009135	PMS2 Sequencing	86009136	PMS2 Deletion, Duplication
Result Code	Result Name							
86009135	PMS2 Sequencing							
86009136	PMS2 Deletion, Duplication							

Chlamydia trachomatis DNA, SDA		
Clinical Significance	C. trachomatis infections are the leading cause of sexually transmitted diseases in the United States. C. trachomatis is known to cause Cervicitis, Pelvic Inflammatory Disease (PID), Infant Conjunctivitis, Infant Pneumonia, Urethritis, Epididymitis and Proctitis. It is also the most frequent cause of non-gonococcal urethritis in men. Among women, the consequences of chlamydial infections are severe if left untreated. Approximately half of chlamydial infections are asymptomatic.	
Effective Date	4/30/2013	
Test Code	17303	
CPT Codes	87491	
Specimen Requirements	20 mL Urine, First void clean catch; Urethral swab or Endocervical swab or Vaginal swab	
Reject Criteria	Only BD Collection Kits are acceptable. Overfilled or under filled Qx UPT; BD ProbeTec Qx CT/GC Amplified Kit w/o swab; BD ProbeTec Qx CT/GC Amp Kit w/cleaning swab. All other specimen types will be rejected.	
Instructions	Urine: Patient should not have urinated for at least 1 hour. Collect the specimen in a sterile, preservative-free collection cup. The patient should collect the first 20-60 mL of voided urine (the first part of the stream - Not Midstream)	

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	<p>into a urine collection cup. Urine should be transferred from collection cup to Qx UPT (Urine Preservative Transport Qx) within 8 hours of collection provided the urine has been stored at 2-30 degrees C. Urine can be held for up to 24 hours prior to transfer to the Qx UPT provided that the urine has been stored at 2-8 degrees C. The correct volume of urine has been added when fluid level is in the fill window. This volume corresponds to 2-3 mL of urine. DO NOT overfill or under fill the tube.</p> <p>Female endocervix: Submit swab in BD ProbeTec CT/GC Qx Amplified Assay Collection Kit for Endocervical specimens.</p> <p>Male urethral: Submit swab in BD ProbeTec CT/GC Qx Amplified Assay Collection Kit for Male Urethral specimens.</p> <p>Vaginal: Submit swab in BD ProbeTec CT/GC Qx Amplified DNA Assay collection kit for Vaginal specimens.</p> <p>If a urine specimen is to be collected, the patient must not have urinated during the previous hour.</p>				
Transport Temperature	See instructions				
Specimen Stability	<p>Urine; Urethral or Endocervical swab Room temperature: 30 Days Refrigerated: 30 Days Frozen: 6 Months</p> <p>Vaginal swab: Room temperature: 14 Days Refrigerated: 14 Days Frozen: 6 Months</p>				
Set-up/Analytic Time	Set-up: Mon-Sat; Report available: 2-3 days				
Reference Range	NOT DETECTED				
Always Message	This test was performed using the BD ProbeTec(TM) Chlamydia trachomatis and Neisseria gonorrhoeae Amplified DNA Assays.				
Methodology	Strand Displacement Amplification				
Assay Category	FDA Approved/Cleared				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>70043703</td> <td>C trachomatis DNA,SDA</td> </tr> </tbody> </table>	Result Code	Result Name	70043703	C trachomatis DNA,SDA
Result Code	Result Name				
70043703	C trachomatis DNA,SDA				

Chlamydia trachomatis DNA, SDA, Pap Vial	
Clinical Significance	C. trachomatis infections are the leading cause of sexually transmitted diseases in the United States. C. trachomatis is known to cause Cervicitis, Pelvic Inflammatory Disease (PID), Infant Conjunctivitis, Infant Pneumonia, Urethritis, Epididymitis and Proctitis. It is also the most frequent cause of non-gonococcal urethritis in men. Among women, the consequences of chlamydial infections are severe if left untreated. Approximately half of chlamydial infections are asymptomatic.
Effective Date	4/30/2013
Test Code	17615
CPT Codes	87491
Specimen Requirements	Liquid Cytology(PreservCyt®) Preservative (ThinPrep®) or TriPath SurePath vials
Reject Criteria	PreservCyt material previously processed for cytology; SurePath material previously processed for cytology; Add-ons to samples for ThinPrep or SurePath; Specimens received after Pap Slide preparation.
Instructions	Liquid Based Cytology: Submit per manufacturer instructions. For samples that have the Pap performed outside of Quest Diagnostics:

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	<p>SurePath(TM) - 0.5 mL pre-aliquot* of SurePath material.** ThinPrep(R) Vial- 0.5 mL pre-aliquot* of PreservCyt(R) material.** *Aliquot before performance of liquid based cytology testing **Refer to manufacturer instructions Transfer aliquot to BD ProbeTec Qx CT/GC LBC diluent tube.</p>					
Transport Temperature	Room temperature					
Specimen Stability	Room temperature: 30 Days Refrigerated: 30 Days Frozen: Unacceptable					
Set-up/Analytic Time	Set-up: Mon-Sat; Report available: 2-3 days					
Reference Range	NOT DETECTED					
Always Message	This test was performed using the BD ProbeTec(TM) Chlamydia trachomatis and Neisseria gonorrhoeae Amplified DNA Assays.					
Methodology	Strand Displacement Amplification					
Assay Category	FDA Approved/Cleared					
Performing Site	Quest Diagnostics Nichols Institute, Valencia					
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>70043705</td> <td>C trachomatis DNA,SDA,Vial</td> </tr> </tbody> </table>		Result Code	Result Name	70043705	C trachomatis DNA,SDA,Vial
Result Code	Result Name					
70043705	C trachomatis DNA,SDA,Vial					

Chlamydia/N. gonorrhoeae DNA, SDA, Pap Vial	
Clinical Significance	<p>C. trachomatis infections are the leading cause of sexually transmitted diseases in the United States. C. trachomatis is known to cause Cervicitis, Pelvic Inflammatory Disease (PID), Infant Conjunctivitis, Infant Pneumonia, Urethritis, Epididymitis and Proctitis. It is also the most frequent cause of non-gonococcal urethritis in men. Among women, the consequences of chlamydial infections are severe if left untreated. Approximately half of chlamydial infections are asymptomatic.</p> <p>Neisseria gonorrhoeae (Gonococci) is the causative agent of gonorrhoeae. In men, this disease generally results in anterior urethritis accompanied by purulent exudate. In women, the disease is most often found in the cervix, but the vagina and uterus may also be affected.</p>
Effective Date	4/30/2013
Test Code	17618
CPT Codes	87491, 87591
Specimen Requirements	Liquid Cytology(PreservCyt®) Preservative (ThinPrep®) or TriPath SurePath vials
Reject Criteria	PreservCyt material previously processed for cytology; SurePath material previously processed for cytology; Add-ons to samples for ThinPrep or SurePath; Specimens received after Pap Slide preparation.
Instructions	<p>Liquid Based Cytology: Submit per manufacturer instructions. For samples that have the Pap performed outside of Quest Diagnostics: SurePath(TM) - 0.5 mL pre-aliquot* of SurePath material.** ThinPrep(R) Vial- 0.5 mL pre-aliquot* of PreservCyt(R) material.** *Aliquot before performance of liquid based cytology testing **Refer to manufacturer instructions Transfer aliquot to BD ProbeTec Qx CT/GC LBC diluent tube.</p>
Transport Temperature	Room temperature
Specimen Stability	Room temperature: 30 Days Refrigerated: 30 Days Frozen: Unacceptable

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Set-up/Analytic Time	Set-up: Mon-Sat; Report available: 2-3 days							
Reference Range	C trachomatis DNA,SDA,Vial: NOT DETECTED N gonorrhoeae DNA,SDA,Vial: NOT DETECTED							
Always Message	This test was performed using the BD ProbeTec(TM) Chlamydia trachomatis and Neisseria gonorrhoeae Amplified DNA Assays.							
Methodology	Strand Displacement Amplification							
Assay Category	FDA Approved/Cleared							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>70043705</td> <td>C trachomatis DNA,SDA,Vial</td> </tr> <tr> <td>70043706</td> <td>N gonorrhoeae DNA,SDA,Vial</td> </tr> </tbody> </table>		Result Code	Result Name	70043705	C trachomatis DNA,SDA,Vial	70043706	N gonorrhoeae DNA,SDA,Vial
Result Code	Result Name							
70043705	C trachomatis DNA,SDA,Vial							
70043706	N gonorrhoeae DNA,SDA,Vial							

Chlamydia/N.gonorrhoeae DNA, SDA	
Clinical Significance	<p>C. trachomatis infections are the leading cause of sexually transmitted diseases in the United States. C. trachomatis is known to cause Cervicitis, Pelvic Inflammatory Disease (PID), Infant Conjunctivitis, Infant Pneumonia, Urethritis, Epididymitis and Proctitis. It is also the most frequent cause of non-gonococcal urethritis in men. Among women, the consequences of chlamydial infections are severe if left untreated. Approximately half of chlamydial infections are asymptomatic.</p> <p>Neisseria gonorrhoeae (Gonococci) is the causative agent of gonorrhoeae. In men, this disease generally results in anterior urethritis accompanied by purulent exudate. In women, the disease is most often found in the cervix, but the vagina and uterus may also be infected.</p>
Effective Date	4/30/2013
Test Code	17305
CPT Codes	87491, 87591
Specimen Requirements	20 mL Urine, First void clean catch; Urethral swab or Endocervical swab or Vaginal swab
Reject Criteria	Only BD Collection Kits are acceptable. Overfilled or under filled Qx UPT; BD ProbeTec Qx CT/GC Amplified Kit w/o swab; BD ProbeTec Qx CT/GC Amp Kit w/cleaning swab. All other specimen types will be rejected.
Instructions	<p>Urine: Patient should not have urinated for at least 1 hour. Collect the specimen in a sterile, preservative-free collection cup.</p> <p>The patient should collect the first 20-60 mL of voided urine (the first part of the stream - Not Midstream) into a urine collection cup. Urine should be transferred from collection cup to Qx UPT (Urine Preservative Transport Qx) within 8 hours of collection provided the urine has been stored at 2-30 degrees C. Urine can be held for up to 24 hours prior to transfer to the Qx UPT provided that the urine has been stored at 2-8 degrees C. The correct volume of urine has been added when fluid level is in the fill window. This volume corresponds to 2-3 mL of urine. DO NOT overfill or under fill the tube.</p> <p>Female endocervix: Submit swab in BD ProbeTec CT/GC Qx Amplified Assay Collection Kit for Endocervical specimens.</p> <p>Male urethral: Submit swab in BD ProbeTec CT/GC Qx Amplified Assay Collection Kit for Male Urethral specimens.</p> <p>Vaginal: Submit swab in BD ProbeTec CT/GC Qx Amplified DNA Assay collection kit for Vaginal specimens.</p> <p>If a urine specimen is to be collected, the patient must not have urinated during the previous hour.</p>
Transport Temperature	See instructions

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Specimen Stability	<p>Urine; Urethral or Endocervical swab Room temperature: 30 Days Refrigerated: 30 Days Frozen: 6 Months</p> <p>Vaginal swab: Room temperature: 14 Days Refrigerated: 14 Days Frozen: 6 Months</p>						
Set-up/Analytic Time	Set-up: Mon-Sat; Report available: 2-3 days						
Reference Range	C trachomatis DNA, SDA: NOT DETECTED N gonorrhoeae DNA, SDA: NOT DETECTED						
Always Message	This test was performed using the BD ProbeTec(TM) Chlamydia trachomatis and Neisseria gonorrhoeae Amplified DNA Assays.						
Methodology	Strand Displacement Amplification						
Assay Category	FDA Approved/Cleared						
Performing Site	Quest Diagnostics Nichols Institute, Valencia						
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>70043703</td> <td>C trachomatis DNA,SDA</td> </tr> <tr> <td>70043704</td> <td>N gonorrhoeae DNA,SDA</td> </tr> </tbody> </table>	Result Code	Result Name	70043703	C trachomatis DNA,SDA	70043704	N gonorrhoeae DNA,SDA
Result Code	Result Name						
70043703	C trachomatis DNA,SDA						
70043704	N gonorrhoeae DNA,SDA						

Neisseria gonorrhoeae DNA, SDA	
Clinical Significance	Neisseria gonorrhoeae (Gonococci) is the causative agent of gonorrhoeae. In men, this disease generally results in anterior urethritis accompanied by purulent exudate. In women, the disease is most often found in the cervix, but the vagina and uterus may also be affected.
Effective Date	4/30/2013
Test Code	17304
CPT Codes	87591
Specimen Requirements	20 mL Urine, First void clean catch; Urethral swab or Endocervical swab or Vaginal swab
Reject Criteria	Only BD Collection Kits are acceptable. Overfilled or under filled Qx UPT; BD ProbeTec Qx CT/GC Amplified Kit w/o swab; BD ProbeTec Qx CT/GC Amp Kit w/cleaning swab. All other specimen types will be rejected.
Instructions	<p>Urine: Patient should not have urinated for at least 1 hour. Collect the specimen in a sterile, preservative-free collection cup. The patient should collect the first 20-60 mL of voided urine (the first part of the stream - Not Midstream) into a urine collection cup. Urine should be transferred from collection cup to Qx UPT (Urine Preservative Transport Qx) within 8 hours of collection provided the urine has been stored at 2-30 degrees C. Urine can be held for up to 24 hours prior to transfer to the Qx UPT provided that the urine has been stored at 2-8 degrees C. The correct volume of urine has been added when fluid level is in the fill window. This volume corresponds to 2-3 mL of urine. DO NOT overfill or under fill the tube.</p> <p>Female endocervix: Submit swab in BD ProbeTec CT/GC Qx Amplified Assay Collection Kit for Endocervical specimens.</p> <p>Male urethral: Submit swab in BD ProbeTec CT/GC Qx Amplified Assay Collection Kit for Male Urethral specimens.</p> <p>Vaginal: Submit swab in BD ProbeTec CT/GC Qx Amplified DNA Assay collection kit for Vaginal specimens.</p> <p>If a urine specimen is to be collected, the patient must not have urinated during the previous hour.</p>

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Transport Temperature	See instructions				
Specimen Stability	<p>Urine; Urethral or Endocervical swab Room temperature: 30 Days Refrigerated: 30 Days Frozen: 6 Months</p> <p>Vaginal swab: Room temperature: 14 Days Refrigerated: 14 Days Frozen: 6 Months</p>				
Set-up/Analytic Time	Set-up: Mon-Sat; Report available: 2-3 days				
Reference Range	NOT DETECTED				
Always Message	This test was performed using the BD ProbeTec(TM) Chlamydia trachomatis and Neisseria gonorrhoeae Amplified DNA Assays.				
Methodology	Strand Displacement Amplification				
Assay Category	FDA Approved/Cleared				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>70043704</td> <td>N gonorrhoeae DNA,SDA</td> </tr> </tbody> </table>	Result Code	Result Name	70043704	N gonorrhoeae DNA,SDA
Result Code	Result Name				
70043704	N gonorrhoeae DNA,SDA				

Neisseria gonorrhoeae DNA, SDA, Pap Vial	
Clinical Significance	Neisseria gonorrhoeae (Gonococci) is the causative agent of gonorrhoeae. In men, this disease generally results in anterior urethritis accompanied by purulent exudate. In women, the disease is most often found in the cervix, but the vagina and uterus may also be affected.
Effective Date	4/30/2013
Test Code	17617
CPT Codes	87591
Specimen Requirements	Liquid Cytology(PreservCyt®) Preservative (ThinPrep®) or TriPath SurePath vials
Reject Criteria	PreservCyt material previously processed for cytology; SurePath material previously processed for cytology; Add-ons to samples for ThinPrep or SurePath; Specimens received after Pap Slide preparation.
Instructions	<p>Liquid Based Cytology: Submit per manufacturer instructions. For samples that have the Pap performed outside of Quest Diagnostics: SurePath(TM) - 0.5 mL pre-aliquot* of SurePath material.** ThinPrep(R) Vial- 0.5 mL pre-aliquot* of PreservCyt(R) material.** *Aliquot before performance of liquid based cytology testing **Refer to manufacturer instructions Transfer aliquot to BD ProbeTec Qx CT/GC LBC diluent tube.</p>
Transport Temperature	Room temperature
Specimen Stability	Room temperature: 30 Days Refrigerated: 30 Days Frozen: Unacceptable
Set-up/Analytic Time	Set-up: Mon-Sat; Report available: 2-3 days
Reference Range	NOT DETECTED
Always Message	This test was performed using the BD ProbeTec(TM) Chlamydia trachomatis and Neisseria gonorrhoeae

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	Amplified DNA Assays.	
Methodology	Strand Displacement Amplification	
Assay Category	FDA Approved/Cleared	
Performing Site	Quest Diagnostics Nichols Institute, Valencia	
CPU Mappings	Result Code	Result Name
	70043706	N gonorrhoeae DNA,SDA,Vial

Test Changes

The following test changes will be effective on the dates indicated below. **Please note that only the information that is changing appears in this update.** Former test names and test codes have been italicized.

Trichomonas vaginalis RNA, Qualitative TMA, Pap Vial		
Effective Date	3/11/2013	
Test Code	90521	
Specimen Requirements	1 mL Liquid Cytology(PreservCyt(R)) Preservative (ThinPrep(R)) in Aptima(R) Vaginal Collection Kit (Orange label) or Aptima(R) Specimen Transfer tube (green label). 0.5 mL (0.5 mL minimum) (SurePath(TM) preservative fluid collected in Aptima(R) Vaginal Collection Kit (Orange label), or Aptima(R) Specimen Transfer tube (green label), or APTIMA(R) Unisex Swab Specimen Collection Kit Original SurePath (TM) vial is no longer acceptable	
Instructions	Labs performing cytology: Aliquot PreservCyt(R)or SurePath(TM) solution before performance of liquid based cytology testing. PreservCyt(R): Transfer 1 mL of PreservCyt(R) solution into APTIMA(R) Vaginal Collection Tube (orange label) or APTIMA(R) Specimen Transfer tube (green label). Ship to lab. SurePath(TM): SurePath fluid must be transferred to APTIMA Specimen Transport Media within 4 days of collection. Tranfer 0.5 mL of SurePath(TM) preservative fluid into APTIMA(R) Vaginal Collection Tube (orange label), APTIMA(R) Specimen Transfer tube (green label), or APTIMA(R) Unisex Swab Specimen Collection Kit. Ship to lab.	
Specimen Stability	Liquid Cytology(PreservCyt(R)) Preservative (ThinPrep(R))	SurePath(TM) preservative fluid
	Room temperature: 14 days Refrigerated: 30 days Frozen: 6 months	Room temperature: 14 days Refrigerated: 14 days Frozen: Not established
Always Message	This test was performed using the APTIMA Trichomonas vaginalis Assay (GEN-PROBE). The performance characteristics of this assay when used to test SurePath(R) specimens have been determined by Quest Diagnostics Nichols Institute, Valencia. Performance characteristics refer to the analytical performance of the test.	
Assay Category	FDA Approved/Cleared/Modified FISH/Molecular assay	
Performing Site	Quest Diagnostics Nichols Institute, Valencia	

Hepatitis C Viral RNA, Quantitative Real-Time PCR w/Rflx to Qualitative TMA	
Effective Date	3/26/2013
Test Code	10051

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Always Message	<p>Please note: the guidelines for the use of new anti- HCV therapies (boceprevir and telaprevir) recommend using a test method that detects plasma viral nucleic acid levels as low as 10 IU/ml. This assay has a lower Limit of Detection of 7.1 IU/ml for genotype 1 and conforms to the recommendation. Quantitation of plasma HCV nucleic acid levels below 43 IU/ml (the Lower Limit of Quantitation for this test) may not be linear, and in this circumstance are reported as "<43 IU/mL HCV RNA Detected".</p> <p>This test was performed using the COBAS(R) AmpliPrep / COBAS(R) TaqMan(R) HCV.</p> <p>http://education.questdiagnostics.com/faq/HCV-RNA-PCR</p> <p>The performance characteristics of this assay have been determined by Quest Diagnostics Nichols Institute, Valencia. Performance characteristics refer to the analytical performance of the test.</p>
Performing Site	Quest Diagnostics Nichols Institute, Valencia

Hepatitis B Virus DNA, Quantitative, Real-Time PCR							
Clinical Significance	Chronic carriers will persist in producing detectable HBV. Patients with chronic liver disease of unknown origin most commonly have HBV that is detected by Viral DNA testing. Quantitative measurement of HBV Viral DNA may be used to monitor progression of disease.						
Effective Date	4/1/2013						
Former Test Name	Hepatitis B DNA Quantitative Real-Time PCR						
Test Code	8137						
Specimen Requirements	3 mL (2.5 mL minimum) Plasma, EDTA (lavender-top) preferred; Plasma, PPT Potassium EDTA (white top) and Serum acceptable.						
Reject Criteria	Unspun PPT tube						
Instructions	Separate plasma from whole blood within 1 day of collection by centrifugation at 800 to 1600 x g for 20 minutes at room temperature. Recommend collecting two tubes per patient to avoid QNS.						
Transport Temperature	Frozen						
Always Message	<p>The method used in this test is Real-Time PCR of the pre-core region of the circular HBV genome.</p> <p>The test was performed using the COBAS(R) AmpliPrep/COBAS(R) TaqMan(R) HBV Test, v2.0 (Roche Molecular Systems, Inc.)</p>						
Performing Site	Quest Diagnostics Nichols Institute, Valencia						
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>210594</td> <td>Hepatitis B Virus DNA</td> </tr> <tr> <td>63203</td> <td>Hepatitis B Virus DNA</td> </tr> </tbody> </table>	Result Code	Result Name	210594	Hepatitis B Virus DNA	63203	Hepatitis B Virus DNA
Result Code	Result Name						
210594	Hepatitis B Virus DNA						
63203	Hepatitis B Virus DNA						
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>2479</td> <td>Hepatitis B Virus MonitR, Chronic</td> </tr> </tbody> </table>	Test Codes:	Name:	2479	Hepatitis B Virus MonitR, Chronic		
Test Codes:	Name:						
2479	Hepatitis B Virus MonitR, Chronic						

Hepatitis C Viral RNA, Quantitative, Real-Time PCR	
Clinical Significance	Quantitative RNA by Real-Time PCR is useful in confirming HCV infection assessing prognosis (prior to the initiation of therapy), prescribing individualized therapy, predicting response to therapy, and monitoring response to therapy.
Effective Date	4/1/2013

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Former Test Name	HCV RNA, Quantitative Real-Time PCR							
Test Code	7577							
Specimen Requirements	3 mL (2.5 mL minimum) Plasma, EDTA (lavender-top) preferred; Plasma, PPT Potassium EDTA (white top) and Serum acceptable.							
Reject Criteria	Received room temperature; Unspun PPT tube; Samples collected using heparin as anticoagulant.							
Instructions	Separate serum or plasma from whole blood within 6 hours of collection by centrifugation at 800-1600 x g for 20 minutes at room temperature. Transfer the plasma to a properly identified, sterile, polypropylene screw-cap vial and ship frozen.							
Transport Temperature	Frozen							
Units Of Measure	HCV RNA, PCR, Quant: IU/mL HCV RNA, PCR, Quant: LogIU/mL							
Always Message	<p>Please note: the guidelines for the use of new anti HCV therapies (boceprevir and telaprevir) recommend using a test method that detects plasma viral nucleic acid levels as low as 10 IU/ml. This assay has a lower Limit of Detection of 7.1 IU/ml for genotype 1 and conforms to the recommendation. Quantitation of plasma HCV nucleic acid levels below 43 IU/ml (the Lower Limit of Quantitation for this test) may not be linear, and in this circumstance are reported as "<43 IU/mL HCV RNA Detected".</p> <p>This test was performed using the COBAS(R) AmpliPrep / COBAS(R) TaqMan(R) HCV Test Kit (Roche Molecular Systems, Inc.).</p> <p>http://education.questdiagnostics.com/faq/HCV-RNA-PCR</p>							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>200357</td> <td>HCV RNA, PCR, Quant</td> </tr> <tr> <td>210944</td> <td>HCV RNA, PCR, Quant</td> </tr> </tbody> </table>		Result Code	Result Name	200357	HCV RNA, PCR, Quant	210944	HCV RNA, PCR, Quant
Result Code	Result Name							
200357	HCV RNA, PCR, Quant							
210944	HCV RNA, PCR, Quant							
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>5906</td> <td>Hepatitis Autoimmune EvaluatR Plus</td> </tr> <tr> <td>RJA</td> <td>Reflex Hepatitis C Viral RNA, Quantitative Real-Time PCR</td> </tr> </tbody> </table>		Test Codes:	Name:	5906	Hepatitis Autoimmune EvaluatR Plus	RJA	Reflex Hepatitis C Viral RNA, Quantitative Real-Time PCR
Test Codes:	Name:							
5906	Hepatitis Autoimmune EvaluatR Plus							
RJA	Reflex Hepatitis C Viral RNA, Quantitative Real-Time PCR							

HIV-1 RNA, Quantitative, Real-Time PCR	
Clinical Significance	<p>This test is intended for use in conjunction with clinical presentation and other laboratory markers of disease progress for the clinical management of HIV-1 infected patients. The test can be used to assess patient prognosis by measuring the baseline HIV-1 RNA level or to monitor the effects of antiretroviral therapy by measuring changes in EDTA plasma HIV-1 RNA levels during the course of antiretroviral treatment.</p> <p>This test is not intended for use as a screening test for the presence of HIV-1 in blood or blood products or as a diagnostic test to confirm the presence of HIV-1 infection.</p>
Effective Date	4/1/2013
Former Test Name	HIV-1 RNA Quantitation [Real Time PCR]
Test Code	7485A
Specimen Requirements	3 mL (2.5 mL minimum) Plasma, EDTA (lavender-top) preferred; Plasma, PPT Potassium EDTA (white top) acceptable.
Reject Criteria	Specimen collected using heparin as anticoagulant; Leaking, uncapped or broken containers; Frozen

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	plasma received in PPT.							
Instructions	Collect plasma in EDTA (lavender-top) or a (white-top) PPT Vacutainer(TM) plasma preparation tube. Separate plasma from the cells by centrifugation within 6 hours after collection. Transfer the plasma to a plastic screw-cap vial and ship frozen.							
Transport Temperature	Frozen							
Always Message	This test was performed using the Cobas(R) AmpliPrep/Cobas(R) Taqman(R) HIV-1 test kit version 2.0 (Roche Molecular Systems, Inc.).							
Performing Site	Quest Diagnostics Nichols Institute, Valencia							
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55132</td> <td>HIV-1 RNA, QN PCR</td> </tr> <tr> <td>55134</td> <td>HIV-1 RNA, QN PCR</td> </tr> </tbody> </table>		Result Code	Result Name	55132	HIV-1 RNA, QN PCR	55134	HIV-1 RNA, QN PCR
Result Code	Result Name							
55132	HIV-1 RNA, QN PCR							
55134	HIV-1 RNA, QN PCR							
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>7485ASR</td> <td>HIV-1 RNA Quantitation [Real Time PCR] w/Serial Reporting</td> </tr> </tbody> </table>		Test Codes:	Name:	7485ASR	HIV-1 RNA Quantitation [Real Time PCR] w/Serial Reporting		
Test Codes:	Name:							
7485ASR	HIV-1 RNA Quantitation [Real Time PCR] w/Serial Reporting							

Sirolimus, LC/MS/MS					
Effective Date	4/1/2013				
Former Test Name	<i>Sirolimus</i>				
Test Code	36712				
Units Of Measure	mcg/L				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>70006</td> <td>Tacrolimus and Sirolimus, LC/MS/MS</td> </tr> </tbody> </table>	Test Codes:	Name:	70006	Tacrolimus and Sirolimus, LC/MS/MS
Test Codes:	Name:				
70006	Tacrolimus and Sirolimus, LC/MS/MS				

Tacrolimus, Highly Sensitive, LC/MS/MS					
Effective Date	4/1/2013				
Test Code	70007				
Reference Range	No definitive therapeutic or toxic ranges have been established. Optimal blood drug levels are influenced by type of transplant, patient response, time post-transplant, co-administration of other drugs, and drug formulation. The following trough range is a suggested guideline: 5.0-20.0 mcg/L				
Units Of Measure	mcg/L				
Performing Site	Quest Diagnostics Nichols Institute, Valencia				
Tests Affected	<table border="1"> <thead> <tr> <th>Test Codes:</th> <th>Name:</th> </tr> </thead> <tbody> <tr> <td>70006</td> <td>Tacrolimus and Sirolimus, LC/MS/MS</td> </tr> </tbody> </table>	Test Codes:	Name:	70006	Tacrolimus and Sirolimus, LC/MS/MS
Test Codes:	Name:				
70006	Tacrolimus and Sirolimus, LC/MS/MS				

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

21-Hydroxylase Antibody [37916X]	
Effective Date	4/8/2013
Test Code	S51345
Reference Range	Negative
Always Message	<p>Patient with indeterminate result may benefit from serial testing to confirm positivity.</p> <p>This test was performed using a kit that has not been approved or cleared by the FDA. The analytical performance characteristics of this test have been determined by Quest Diagnostics Nichols Institute, San Juan Capistrano. This test should not be used for diagnosis without confirmation by other medically established means.</p>
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano

CBC (Includes Differential and Platelets)																															
Effective Date	4/8/2013																														
Test Code	6399																														
Specimen Stability	Refrigerated: 48 hours																														
Reference Range	<p>Remove the reference ranges from the following components within the CBC:</p> <ul style="list-style-type: none"> % Neutrophilic Bands % Neutrophils % Lymphocytes % Monocytes % Eosinophils % Basophils <p>Reference ranges will continue to be reported for the absolute cell counts.</p>																														
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano																														
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>30000000</td> <td>WBC</td> </tr> <tr> <td>30000100</td> <td>RBC</td> </tr> <tr> <td>30000200</td> <td>HGB</td> </tr> <tr> <td>30000300</td> <td>HCT</td> </tr> <tr> <td>30000400</td> <td>MCV</td> </tr> <tr> <td>30000500</td> <td>MCH</td> </tr> <tr> <td>30000600</td> <td>MCHC</td> </tr> <tr> <td>30000800</td> <td>PLT</td> </tr> <tr> <td>30004600</td> <td>MPV</td> </tr> <tr> <td>30000700</td> <td>RDW</td> </tr> <tr> <td>30001700</td> <td>Absolute Neutrophils</td> </tr> <tr> <td>30001110</td> <td>Absolute Band Neutrophils</td> </tr> <tr> <td>30001310</td> <td>Absolute Metamyelocytes</td> </tr> <tr> <td>30001510</td> <td>Absolute Myelocytes</td> </tr> </tbody> </table>	Result Code	Result Name	30000000	WBC	30000100	RBC	30000200	HGB	30000300	HCT	30000400	MCV	30000500	MCH	30000600	MCHC	30000800	PLT	30004600	MPV	30000700	RDW	30001700	Absolute Neutrophils	30001110	Absolute Band Neutrophils	30001310	Absolute Metamyelocytes	30001510	Absolute Myelocytes
Result Code	Result Name																														
30000000	WBC																														
30000100	RBC																														
30000200	HGB																														
30000300	HCT																														
30000400	MCV																														
30000500	MCH																														
30000600	MCHC																														
30000800	PLT																														
30004600	MPV																														
30000700	RDW																														
30001700	Absolute Neutrophils																														
30001110	Absolute Band Neutrophils																														
30001310	Absolute Metamyelocytes																														
30001510	Absolute Myelocytes																														

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	30001530	Absolute Promyelocytes
	30002110	Absolute Lymphocytes
	30002400	Absolute Monocytes
	30002700	Absolute Eosinophils
	30003000	Absolute Basophils
	30003500	Absolute Blasts
	30003610	Absolute Nucleated RBC
	30000900	Neutrophils
	30001100	Band Neutrophils
	30001300	Metamyelocytes
	30001500	Myelocytes
	30001520	Promyelocytes
	30001800	Lymphocytes
	30002000	Reactive Lymphocytes
	30002200	Monocytes
	30002500	Eosinophils
	30002800	Basophils
	30003400	Blasts
	30003600	Nucleated RBC
	30004200	Comment(s)
Additional Information	Update report format.	

Dihydrotestosterone, LC/MS/MS [90567]	
Effective Date	4/8/2013
Test Code	S51491
Specimen Requirements	1 mL (0.6 mL minimum) serum collected in a red-top tube (no gel)
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano

Estriol, LC/MS/MS, Serum [34883]	
Effective Date	4/8/2013
Former Test Name	Estriol Serum Total [22517P]
Test Code	S49575
Specimen Requirements	0.5 mL (0.3 mL minimum) serum collected in red-top tube (no gel) Plasma is unacceptable

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano	
Tests Affected	Test Codes:	Name:
	36742	Estrogens, Fractionated, LC/MS/MS

Muramidase (Lysozyme), Serum		
Effective Date	4/8/2013	
Former Test Name	Lysozyme, Serum [619X]	
Test Code	S51290	
Specimen Requirements	2 mL (0.5 minimum) serum	
Instructions	Centrifuge serum specimens within 1 hour of collection. Transfer serum to sterile, plastic, screw-capped vial(s). Freeze and ship frozen samples on dry ice.	
Transport Temperature	Frozen	
Set-up/Analytic Time	Set up: Mon, Thurs; Report available: 1-5 days	
Reference Range	7.0-15.0 mcg/mL	
Methodology	Enzymatic	
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano	
CPU Mappings	Result Code	Result Name
	110202	Muramidase (Lysozyme), Serum

Porphobilinogen, Quantitative, 24-Hour Urine			
Effective Date	4/8/2013		
Former Test Name	Porphobilinogen 24hr Urine [7294N]		
Test Code	S51657		
Always Message	Remove Always Message		
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano		
CPU Mappings	Result Code	Type	Result Name
	110939	AOE	Total Volume
	110940		Porphobilinogen, 24 hr Ur
	86008983		Interpretation
Additional Information	Update report format to add interpretation analyte.		
Tests Affected	Test Codes:	Name:	
	S48559	Porphyrins, Fractionated, Quantitative and Porphobilinogen, 24-Hour Urine	

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Porphobilinogen, Quantitative, Random Urine							
Effective Date	4/8/2013						
Former Test Name	Porphobilinogen Random Urine [84889N]						
Test Code	S51658						
Always Message	Remove Always Message						
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano						
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>110941</td> <td>Porphobilinogen, Rand Ur</td> </tr> <tr> <td>86008984</td> <td>Interpretation</td> </tr> </tbody> </table>	Result Code	Result Name	110941	Porphobilinogen, Rand Ur	86008984	Interpretation
	Result Code	Result Name					
	110941	Porphobilinogen, Rand Ur					
86008984	Interpretation						
Additional Information	Update report format to add interpretation analyte.						

Porphyrins, Fractionated, Quantitative, Random Urine	
Effective Date	4/8/2013
Former Test Name	Porphyrins Fractionated Urine Random [9324N]
Test Code	S48558
Specimen Requirements	2 mL (1.5 mL minimum) random urine, no preservative
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano

Porphyrins, Total, Plasma							
Effective Date	4/8/2013						
Former Test Name	Porphyrins, Total Plasma [10290X]						
Test Code	S51406						
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano						
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>110328</td> <td>Total Porphyrins</td> </tr> <tr> <td>86009304</td> <td>Interpretation</td> </tr> </tbody> </table>	Result Code	Result Name	110328	Total Porphyrins	86009304	Interpretation
	Result Code	Result Name					
	110328	Total Porphyrins					
86009304	Interpretation						
Additional Information	Update report format to add interpretation analyte.						

Protein C Antigen	
Effective Date	4/8/2013
Test Code	4948
Reject Criteria	Received room temperature; received refrigerated; gross hemolysis
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Tests Affected	Test Codes:	Name:
	8757	Protein C Activity and Antigen

Thyroid Abs Evaluation											
Effective Date	4/16/2013										
Test Code	3060										
Always Message	<p>Remove: Reference Range: <140% of baseline of Reference Control</p> <p>Thyroid Stimulating Immunoglobulins (TSI) can engage the TSH receptors resulting in hyperthyroidism in Graves' disease patients. TSI levels can be useful in monitoring the clinical outcome of Graves' disease as well as assessing the potential for hyperthyroidism from maternal-fetal transfer. TSI results greater than or equal to (>=) 140% of the Reference Control are considered positive.</p> <p>NOTE: A serum TSH level greater than 350 micro-International Units/mL can interfere with the TSI bioassay and potentially give false positive results.</p> <p>Patients who are pregnant and are suspected of having hyperthyroidism should have both a TSI and human Chorionic Gonadotropin (hCG) tests measured. A serum hCG level greater than 40,625 mIU/mL can interfere with the TSI bioassay and may give false negative results. In these patients it is recommended that a second TSI is obtained when the hCG concentration falls below 40,625 mIU/mL (usually after approximately 20-weeks gestation).</p>										
Performing Site	Quest Diagnostics Nichols Institute, Valencia										
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>11165</td> <td>Thyroglobulin Autoantibodies</td> </tr> <tr> <td>30876</td> <td>Thyroid Peroxidase Autoantibodies</td> </tr> <tr> <td>32490</td> <td>TBII</td> </tr> <tr> <td>32525</td> <td>TSH</td> </tr> </tbody> </table>	Result Code	Result Name	11165	Thyroglobulin Autoantibodies	30876	Thyroid Peroxidase Autoantibodies	32490	TBII	32525	TSH
Result Code	Result Name										
11165	Thyroglobulin Autoantibodies										
30876	Thyroid Peroxidase Autoantibodies										
32490	TBII										
32525	TSH										
Additional Information	Remove TSI from test offering, result code 32516.										

Varicella-Zoster Virus Antibody (IgG)	
Effective Date	4/30/2013
Test Code	8761
Reject Criteria	Gross hemolysis; gross lipemia; grossly icteric
Reference Range	<p>< or = 0.90</p> <p>Interpretive Criteria < or = 0.90 Negative-No VZV IgG Antibody Detected 0.91-1.09 Equivocal > or = 1.10 Positive-VZV IgG Antibody Detected</p>
Always Message	A positive result indicates that the patient has antibody to VZV but does not differentiate between an active or past infection. The clinical diagnosis must be interpreted in conjunction with the clinical signs and symptoms of the patient. This assay reliably measures immunity due to previous infection but may not be sensitive enough to detect antibodies induced by vaccination. Thus, a negative result in a vaccinated individual does not necessarily indicate susceptibility to VZV infection.

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

Performing Site	Quest Diagnostics Nichols Institute, Valencia	
Tests Affected	Test Codes:	Name:
	8756	Varicella-Zoster Virus Antibodies (IgG, IgM)

Redirects

TSI (Thyroid Stimulating Immunoglobulin)		
Effective Date	4/16/2013	
Former Test Code	1092	
Test Code	30551	
Reject Criteria	Gross hemolysis; Gross lipemia; Gross icterea	
Set-up/Analytic Time	Set-up: Tues-Sat; Report available: 1-5 days	
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano	
CPU Mappings	Result Code	Result Name
	85986240	TSI

Vitamin B1 (Thiamine), Blood, LC/MS/MS	
Clinical Significance	Vitamin B1 is required for branched-chain amino acid and carbohydrate metabolism. Vitamin B1 deficiency is most often due to alcoholism or chronic illness. In the early stage, patients with Vitamin B1 deficiency exhibit anorexia, irritability, apathy, and generalized weakness. Prolonged deficiency causes beriberi.
Effective Date	4/23/2013
Former Test Name	Vitamin B1 (Thiamine) Whole Blood
Former Test Code	3515W
Test Code	5042X
Specimen Requirements	2 mL (1 mL) whole blood collected in an EDTA (lavender-top) tube
Reject Criteria	Received room temperature; Received refrigerated
Instructions	Transfer whole blood to a plastic shipping vial to prevent breakage. Wrap tube in aluminum foil to protect from light. Freeze immediately. Samples received frozen but not protected from light are acceptable. Wrap sample in foil before sending to lab, no disclaimer is required. Samples received refrigerated but not protected from light are unacceptable and should be rejected. Overnight fasting is recommended
Specimen Stability	Room temperature: 72 hours Refrigerated: 5 days Frozen: 5 weeks
Set-up/Analytic Time	Set up: Mon-Fri; Report available: 3-6 days

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	This assay has an overnight incubation period					
Reference Range	78-185 nmol/L					
Units Of Measure	nmol/L					
Methodology	Liquid Chromatography/Tandem Mass Spectrometry					
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano					
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>85987120</td> <td>Vitamin B1, Blood</td> </tr> </tbody> </table>		Result Code	Result Name	85987120	Vitamin B1, Blood
Result Code	Result Name					
85987120	Vitamin B1, Blood					

HE4, Ovarian Cancer Monitoring					
Effective Date	5/6/2013				
Former Test Code	2040				
Test Code	16500				
Specimen Requirements	0.5 mL (0.3 mL minimum) serum				
Reject Criteria	Gross hemolysis; lipemia; icteria; presence of fibrin, red blood cells, or other particulate matter; obvious microbial contamination				
Instructions	Sample collected in red-top tube must be allowed to clot for at least 30 minutes but no longer than 1 hour.				
Transport Temperature	Frozen				
Specimen Stability	Room temperature: 48 hours Refrigerated: 72 hours Frozen: 28 days				
Set-up/Analytic Time	Set up: Mon, Wed, Fri; Report available: 2-5 days				
Reference Range	Female Premenopausal < or = 70 pmol/L Female Postmenopausal < or = 140 pmol/L				
Units Of Measure	pmol/L				
Methodology	Chemiluminescent Microparticle Immunoassay (CMIA)				
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano.				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>55005776</td> <td>HE4, Ovarian CA Monitoring</td> </tr> </tbody> </table>	Result Code	Result Name	55005776	HE4, Ovarian CA Monitoring
Result Code	Result Name				
55005776	HE4, Ovarian CA Monitoring				

Discontinued

Hepatitis C Antibody w/Reflex RIBA	
Effective Date	4/16/2013
Test Code	2445
Additional Information	The recommended alternative is 7493 Hepatitis C Antibody w/Rfx to RNA Quantitation performed at Quest Diagnostics

March 2013 - Monthly Update, Quest Diagnostics Nichols Institute, Valencia

	Nichols Institute, Valencia	
Tests Affected	Test Codes:	Name:
	2445B	Hepatitis C Antibody w/Reflex RIBA + Bands

Thyroid Stimulating Immunoglobulins w/TSH	
Effective Date	4/16/2013
Test Code	1091
Performing Site	Quest Diagnostics Nichols Institute, Valencia and San Juan Capistrano
Additional Information	The recommended alternatives are: 3250 TSH performed at Quest Diagnostics Nichols Institute, Valencia and 30551 TSI (Thyroid Stimulating Immunoglobulin) performed at Quest Diagnostics Nichols Institute, San Juan Capistrano

Chlamydia/N. gonorrhoeae DNA, SDA	
Effective Date	4/30/2013
Test Code	S51330
Additional Information	The replacement test is 17305 Chlamydia/N. gonorrhoeae DNA, SDA performed at Quest Diagnostics Nichols Institute, Valencia

New York Update

Herpes Simplex Virus, Type 1/2 DNA, Real-Time PCR, PAP Vial	
Message	Test is now approved for New York patients.
Effective Date	3/26/2013
Test Code	90569
Performing Site	Quest Diagnostics Nichols Institute, Valencia