

10/28/2015 - Immediate Action, Quest Diagnostics Nichols Institute, Valencia

TEST CHANGES				
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.				
Test Code	Former Test Code	Test Name	Effective Date	Page #
<u>92777</u>		QNatal™ Advanced	10/14/2015	1
<u>92688</u>		Liver Fibrosis, FibroTest-ActiTest Panel	11/9/2015	1

SEND OUTS				
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.				
Test Code	Former Test Code	Test Name	Effective Date	Page #
S51264		Fatty Acid Profile, Essential (C12-C22), Serum	10/28/2015	2
<u>91810</u>		Rabies RFFIT Serology End Point Titer	10/28/2015	3
S52235		Rabies Vaccine Response End-Point Titer	10/28/2015	3
<u>11194</u>		14-3-3 Protein, Spinal Fluid	11/5/2015	3

Due to the background information related to these changes, we are unable to use our normal method of communication. These changes listed in this document are effective in less than 30 days. Please note the individual effective dates below, as some of these changes require IMMEDIATE ACTION.

Test Changes

The following test changes will be effective on the dates indicated below. Please note information that is changing appears in bold text in this update. Former test names and test codes have been italicized.

QNatal™ Advanced	
Message	**This test is not available for California or New York patient testing**
Effective Date	10/14/2015
Test Code	<i>92777</i>
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano

Liver Fibrosis, FibroTest-ActiTest Panel																							
Effective Date	11/9/2015																						
Test Code	<i>92688</i>																						
Reference Range	<table border="1"> <thead> <tr> <th>Result Name</th> <th>Reference Range</th> </tr> </thead> <tbody> <tr> <td>Fibrosis Score</td> <td>No reference range available</td> </tr> <tr> <td>Fibrosis Stage</td> <td>No reference range available</td> </tr> <tr> <td>Fibrosis Interpretation</td> <td> <table border="0"> <tr> <td>FibroTest Score</td> <td>Metavir Score</td> </tr> <tr> <td>0.00-0.21</td> <td>F0 no fibrosis</td> </tr> <tr> <td>0.22-0.27</td> <td>F0-F1</td> </tr> <tr> <td>0.28-0.31</td> <td>F1 minimal fibrosis</td> </tr> <tr> <td>0.32-0.48</td> <td>F1-F2</td> </tr> <tr> <td>0.49-0.58</td> <td>F2 moderate fibrosis</td> </tr> <tr> <td>0.59-0.72</td> <td>F3 advanced fibrosis</td> </tr> </table> </td> </tr> </tbody> </table>	Result Name	Reference Range	Fibrosis Score	No reference range available	Fibrosis Stage	No reference range available	Fibrosis Interpretation	<table border="0"> <tr> <td>FibroTest Score</td> <td>Metavir Score</td> </tr> <tr> <td>0.00-0.21</td> <td>F0 no fibrosis</td> </tr> <tr> <td>0.22-0.27</td> <td>F0-F1</td> </tr> <tr> <td>0.28-0.31</td> <td>F1 minimal fibrosis</td> </tr> <tr> <td>0.32-0.48</td> <td>F1-F2</td> </tr> <tr> <td>0.49-0.58</td> <td>F2 moderate fibrosis</td> </tr> <tr> <td>0.59-0.72</td> <td>F3 advanced fibrosis</td> </tr> </table>	FibroTest Score	Metavir Score	0.00-0.21	F0 no fibrosis	0.22-0.27	F0-F1	0.28-0.31	F1 minimal fibrosis	0.32-0.48	F1-F2	0.49-0.58	F2 moderate fibrosis	0.59-0.72	F3 advanced fibrosis
Result Name	Reference Range																						
Fibrosis Score	No reference range available																						
Fibrosis Stage	No reference range available																						
Fibrosis Interpretation	<table border="0"> <tr> <td>FibroTest Score</td> <td>Metavir Score</td> </tr> <tr> <td>0.00-0.21</td> <td>F0 no fibrosis</td> </tr> <tr> <td>0.22-0.27</td> <td>F0-F1</td> </tr> <tr> <td>0.28-0.31</td> <td>F1 minimal fibrosis</td> </tr> <tr> <td>0.32-0.48</td> <td>F1-F2</td> </tr> <tr> <td>0.49-0.58</td> <td>F2 moderate fibrosis</td> </tr> <tr> <td>0.59-0.72</td> <td>F3 advanced fibrosis</td> </tr> </table>	FibroTest Score	Metavir Score	0.00-0.21	F0 no fibrosis	0.22-0.27	F0-F1	0.28-0.31	F1 minimal fibrosis	0.32-0.48	F1-F2	0.49-0.58	F2 moderate fibrosis	0.59-0.72	F3 advanced fibrosis								
FibroTest Score	Metavir Score																						
0.00-0.21	F0 no fibrosis																						
0.22-0.27	F0-F1																						
0.28-0.31	F1 minimal fibrosis																						
0.32-0.48	F1-F2																						
0.49-0.58	F2 moderate fibrosis																						
0.59-0.72	F3 advanced fibrosis																						

	0.73-0.74 0.75-1.00	F3-F4 F4 severe fibrosis
Necroinflammat Act Score	No reference range available	
Necroinflammat Act Grade	No reference range available	
Necroinflammat Interp	ActiTest Score	Metavir Score
	0.00-0.17	A0 no activity
	0.18-0.29	A0-A1
	0.30-0.36	A1 minimal activity
	0.37-0.52	A1-A2
	0.53-0.60	A2 significant activity
	0.61-0.62	A2-A3
	0.63-1.00	A3 severe activity
Alpha-2-Macroglobulin - No change		
Haptoglobin- No change		
Apolipoprotein A1 - No change		
Total Bilirubin - No change		
GGT - No change		
ALT - No change		
Performing Site	Quest Diagnostics Nichols Institute, San Juan Capistrano	

Test Send Out (Referrals)

Fatty Acid Profile, Essential (C12-C22), Serum			
Effective Date	10/28/2015		
Test Code	S51264		
Specimen Requirements	0.5 mL (0.15 mL minimum) serum		
Interface Mapping	Result Code	Result Name	Unit Of Measure
	109958	Lauric Acid, C12:0	nmol/mL
	109959	Myristic Acid, C14:0	nmol/mL
	109960	Hexadecenoic Acid, C16:1w9	nmol/mL
	109961	Palmitoleic Acid, C16:1w7	nmol/mL
	109962	Palmitic Acid, C16:0	nmol/mL
	109963	g-Linolenic Acid, C18:3w6	nmol/mL
	109964	a-Linolenic Acid, C18:3w3	nmol/mL
	109965	Linoleic Acid, C18:2w6	nmol/mL
	109966	Oleic Acid, C18:1w9	nmol/mL
	109967	Vaccenic Acid, C18:1w7	nmol/mL
	109968	Stearic Acid, C18:0	nmol/mL
	109969	EPA, C20:5w3	nmol/mL
	109970	Arachidonic Acid, C20:4w6	nmol/mL
	109971	Mead Acid, C20:3w9	nmol/mL
	109984	h-g-Linolenic C20:3w6	nmol/mL
	109973	Arachidic Acid, C20:0	nmol/mL

109974	DHA, C22:6w3	nmol/mL
109975	DPA, C22:5w6	nmol/mL
109976	DPA, C22:5w3	nmol/mL
109977	DTA, C22:4w6	nmol/mL
109978	Docosenoic Acid, C22:1	nmol/mL
109993	Nervonic Acid, C24:1w9	nmol/mL
109980	Triene Tetraene Ratio	
109981	Total Saturated Acid	mmol/L
109982	Total Monounsaturated Acid	mmol/L
109983	Total Polyunsaturated Acid	mmol/L
109985	Total w3	mmol/L
109986	Total w6	mmol/L
109987	Total Fatty Acids	mmol/L
109988	Interpretation	

Rabies RFFIT Serology End Point Titer							
Effective Date	10/28/2015						
Test Code	91810						
Specimen Stability	<table border="1"> <tr> <td>Room temperature</td> <td>48 hours</td> </tr> <tr> <td>Refrigerated</td> <td>14 days</td> </tr> <tr> <td>Frozen</td> <td>28 days</td> </tr> </table>	Room temperature	48 hours	Refrigerated	14 days	Frozen	28 days
Room temperature	48 hours						
Refrigerated	14 days						
Frozen	28 days						
Additional Information	This test is being reactivated at Quest Diagnostics, Nichols Institute, Valencia.						

Rabies Vaccine Response End-Point Titer	
Message	**This test is not available for Florida patient testing. Please use test code 91810- Rabies RFFIT Serology End Point Titer for Florida patient testing. **
Effective Date	10/28/2015
Test Code	S52235

14-3-3 Protein, Spinal Fluid					
Effective Date	11/5/2015				
Test Code	11194				
Reference Range	<table border="1"> <thead> <tr> <th>Result Name</th> <th>Reference Range</th> </tr> </thead> <tbody> <tr> <td>14.3.3 Protein, CSF</td> <td>Normal: < or = 2.0 ng/mL Elevated: > 2.0 ng/mL</td> </tr> </tbody> </table>	Result Name	Reference Range	14.3.3 Protein, CSF	Normal: < or = 2.0 ng/mL Elevated: > 2.0 ng/mL
Result Name	Reference Range				
14.3.3 Protein, CSF	Normal: < or = 2.0 ng/mL Elevated: > 2.0 ng/mL				