

8/12/2013 - New Release, Quest Diagnostics Nichols Institute, Valencia

NEW TESTS			
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.			
Test Code	Test Name	Effective Date	Page #
<u>91830</u>	Acetyl Fentanyl, Qualitative, Blood	8/19/2013	1
<u>91831</u>	Acetyl Fentanyl, Qualitative, Serum or Plasma	8/19/2013	1
<u>91829</u>	Acetyl Fentanyl, Qualitative, Urine	8/19/2013	2

New Test Offerings

The following tests will be available through Quest Diagnostics on the dates indicated below.

Acetyl Fentanyl, Qualitative, Blood					
Clinical Significance	Acetyl fentanyl is a fentanyl analog previously undocumented in illicit drug use that is up to five times more potent than heroin. Acetyl fentanyl's enhanced potency has led to a number of reported deaths due to overdose. The Centers for Disease Control (CDC) recommends laboratories to provide analytical options for acetyl fentanyl in order to provide enhanced surveillance for healthcare providers.				
Effective Date	8/19/2013				
Test Code	91830				
CPT Codes	80101				
Specimen Requirements	<p>Preferred: 3 mL (1 mL minimum) whole blood collected in a sodium fluoride (gray-top) tube</p> <p>Acceptable: 3 mL (1 mL minimum) whole blood collected in a sodium heparin (green-top) tube or EDTA (lavender-top) tube</p>				
Transport Temperature	Room temperature				
Specimen Stability	<p>Room temperature: 5 days</p> <p>Refrigerated: 7 days</p> <p>Frozen: 30 days</p>				
Set-up/Analytic Time	Set up: Mon-Fri; Report available: 2 days				
Reference Range	Negative				
Always Message	Reference Range: Negative				
Methodology	Liquid Chromatography/Tandem Mass Spectrometry (QTrap)				
Assay Category	Laboratory Developed Test				
Performing Site	Quest Diagnostics Nichols Institute, Chantilly				
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86010026</td> <td>Acetyl Fentanyl, QL, Blood</td> </tr> </tbody> </table>	Result Code	Result Name	86010026	Acetyl Fentanyl, QL, Blood
Result Code	Result Name				
86010026	Acetyl Fentanyl, QL, Blood				

Acetyl Fentanyl, Qualitative, Serum or Plasma	
Clinical Significance	Acetyl fentanyl is a fentanyl analog previously undocumented in illicit drug use that is up to five times more potent than

8/12/2013 - New Release, Quest Diagnostics Nichols Institute, Valencia

	heroin. Acetyl fentanyl's enhanced potency has led to a number of reported deaths due to overdose. The Centers for Disease Control (CDC) recommends laboratories to provide analytical options for acetyl fentanyl in order to provide enhanced surveillance for healthcare providers.					
Effective Date	8/19/2013					
Test Code	91831					
CPT Codes	80101					
Specimen Requirements	Preferred: 3 mL (1 mL minimum) serum collected in a red-top tube (no gel) Acceptable: 3 mL plasma (1 mL minimum) collected in an EDTA (lavender-top) tube					
Reject Criteria	Serum separator tube					
Transport Temperature	Room temperature					
Specimen Stability	Room temperature: 5 days Refrigerated: 7 days Frozen: 30 days					
Set-up/Analytic Time	Set up: Mon-Fri; Report available: 2 days					
Reference Range	Negative					
Always Message	Reference Range: Negative					
Methodology	Liquid Chromatography/Tandem Mass Spectrometry (QTrap)					
Assay Category	Laboratory Developed Test					
Performing Site	Quest Diagnostics Nichols Institute, Chantilly					
CPU Mappings	<table border="1"> <thead> <tr> <th>Result Code</th> <th>Result Name</th> </tr> </thead> <tbody> <tr> <td>86010027</td> <td>Acetyl Fentanyl, QL, S/P</td> </tr> </tbody> </table>		Result Code	Result Name	86010027	Acetyl Fentanyl, QL, S/P
Result Code	Result Name					
86010027	Acetyl Fentanyl, QL, S/P					

Acetyl Fentanyl, Qualitative, Urine	
Clinical Significance	Acetyl fentanyl is a fentanyl analog previously undocumented in illicit drug use that is up to five times more potent than heroin. Acetyl fentanyl's enhanced potency has led to a number of reported deaths due to overdose. The Centers for Disease Control (CDC) recommends laboratories to provide analytical options for acetyl fentanyl in order to provide enhanced surveillance for healthcare providers.
Effective Date	8/19/2013
Test Code	91829
CPT Codes	80101
Specimen Requirements	5 mL (1 mL minimum) random urine submitted in a leakproof urine container
Transport Temperature	Room temperature
Specimen Stability	Room temperature: 5 days Refrigerated: 7 days Frozen: 30 days
Set-up/Analytic Time	Set up: Mon-Fri; Report available: 2 days
Reference Range	Negative
Always Message	Reference Range: Negative

Methodology	Liquid Chromatography/Tandem Mass Spectrometry (QTrap)	
Assay Category	Laboratory Developed Test	
Performing Site	Quest Diagnostics Nichols Institute, Chantilly	
CPU Mappings	Result Code	Result Name
	86010025	Acetyl Fentanyl, QL, Urine