

12/2/2014 - TSO Update, Quest Diagnostics Nichols Institute, Valencia

TEST CHANGES				
Please Note: Not all test codes assigned to each assay are listed in the table of contents. Please refer to the complete listing on the page numbers indicated.				
Test Code	Former Test Code	Test Name	Effective Date	Page #
<u>S47310</u>		Acylglycines Quantitative, Urine	12/8/2014	1

Due to these assays being performed by outside vendors, we are unable to use our normal method of communication. Some of the changes listed in this document may be effective in less than 30 days. Please note the individual effective dates below, as these changes may require **IMMEDIATE ACTION**.

## Test Changes

Acylglycines Quantitative, Urine	
Effective Date	12/8/2014
Test Code	S47310
Specimen Requirements	Preferred: 10 mL <b>random urine</b> collected in a sterile screw cap container  <b>Remove: 24-hour urine as an acceptable specimen</b>
Specimen Stability	Room temperature: Unacceptable <b>Refrigerated: 9 days</b> <b>Frozen: 416 days</b>
Set-up/Analytic Time	Set up: Mon, Wed, Fri; <b>Report available: 5-14 days</b>