

Specialty Immunohistochemistry(IHC) Update:

The following IHC stains have been added to our menu. For a complete listing of all IHC testing please go to www.specialtylabs.com or contact Client Relations at 800-421-4449.

Specimen/Stability	Formalin fixed, paraffin-embedded tissue; Ambient
Collection Notes	Shipping on cold pack is recommended during the warm weather months. Ship by overnight courier to <i>Specialty</i> .
Clinical Utility	Useful for identification of antigens present in paraffin-embedded tissues. Properly selected panels may aid in the identification of tumor type and subclassification. Additional markers may be useful in assessing proliferation (Ki67) or prognosis (e.r., bcl-2 or HER-2/ <i>neu</i>) in selected tumor types.
Schedule	Monday-Sunday
Reported	1-3 days
Notes	Diagnostic code required for third party reimbursement. Please send copy of pathology report and other applicable test results. <i>Specialty</i> continues to enhance client options for IHC staining; please contact our Oncology/IHC department at 888-567-5225 regarding specific IHC requests not listed on this menu.

IHC205	WT-1	IHC224	Calponin B	Retic
IHC218	CD23	IHC222	Placental	PAS
IHC216	CD35		Lactogen	PAS with amylase
IHC207	MOC31	IHC225	TIA-1	
IHC214	p504s	IHC232	Ham56	<u>Coming in August</u>
IHC211	Cathespin D	IHC228	A1-ACT	PMS2
IHC215	Glycophorin A	IHC233	A-AT	MGMT
IHC217	CEA	IHC231	HMB-45 RED	MDM2
	polyclonal	IHC234	Melan A RED	
IHC213	Pgp9.5	IHC227	PIN 4 Cocktail	
IHC219	Collagen Iv		Double Stain	
IHC220	Mast Cell			
	Tryptase	<u>Coming in June</u>		
IHC206	CD7	CD25		
IHC208	IgA	<i>H.Pylori</i>		
IHC209	IgG	HBcAg		
IHC210	IgM	Herpes (HSV)		
IHC226	MNDA	MLH-1		
IHC221	MUM 1	MHS-2		
IHC212	pS2	MSH-6		
IHC229	CD61	<i>P. carinii</i>		
IHC230	RCC	P27		
IHC235	hepatocyte	Kappa/Lambda double stain		
IHC223	Cyclin E	Fascin		